

# Assay Training: Xpert<sup>®</sup> Carba-R

Cepheid Technical Support  
CE-IVD



 **Cepheid<sup>®</sup>**  
*A better way.*

# Training Agenda

- **Xpert Carba-R Training**
  - Reagents
  - Sample collection
  - Kit storage and handling
  - Precautions
  - Preparing cartridge
- **Quality Control**
- **Results analysis**
- **Discussion and Q&A**



# Xpert Carba-R Training Objectives


**At the end of the training, user will be able to:**

- Properly store and handle the Xpert Carba-R cartridge kits.
- Follow proper laboratory safety precautions.
- Collect appropriate specimen types and transport specimen.
- Perform the cartridge set up and run the assay.
- Report the various software-generated results.
- Understand assay control strategy.

# Introduction on disease state

## *What is Carbapenem resistance (Carba-R)?*

- **Carbapenems are a type of beta-lactam antibiotic, often an antibiotic of last resort.**
  - **Carbapenem resistance (Carba-R) is a result of chromosomal and highly mobile plasmid mediated resistance genes associated with gram-negative enteric organisms, Pseudomonas, and Acinetobacters.**
  - **Infection with carbapenem non-susceptible organisms are associated with high mortality (up to 40-50%).<sup>1</sup>**
  - **Screening for colonized patients can expedite infection control practices and reduce likelihood of nosocomial spread.**



<sup>1</sup>Guidance for control of carbapenem-resistant Enterobacteriaceae (CRE): 2012 CRE Toolkit. Division of Healthcare Quality Promotion. Centers for Disease Control and Prevention. Atlanta, Georgia. 2012.

# The Cepheid Solution



- **Simultaneous detection and differentiation**
  - Five Carba-R gene targets
  - Two internal controls for each individual sample
  - Sample Processing Control (SPC)
  - Probe Check Control (PCC)
- **High sensitivity and specificity**
- **Simple and easy to use**
  - Closed cartridge system
- **Results in 48 minutes**
- **On-demand results 24/7**
- **Random access**

# Intended Use

The Cepheid Xpert Carba-R Assay, performed on the GeneXpert® Instrument Systems, is a qualitative *in vitro* diagnostic test designed for rapid detection and differentiation of the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP-1</sub> gene sequences associated with carbapenem-non-susceptibility in gram-negative bacteria obtained from rectal swab specimens from patients at risk for intestinal colonization with carbapenem-non-susceptible bacteria. The test utilizes automated real-time polymerase chain reaction (PCR). The Xpert Carba-R Assay is intended to aid in the detection of carbapenem-non susceptible bacteria that colonize patients in healthcare settings. The Xpert Carba-R Assay is not intended to guide or monitor treatment for carbapenem-non-susceptible bacterial infections.

Concomitant cultures are necessary to recover organisms for epidemiological typing, antimicrobial susceptibility testing, and for further confirmatory identification of carbapenem-non-susceptible bacteria.

# System and Reagent Requirements

## GeneXpert Systems

- **6 color modules**
- **GX DX Software v4.3 or higher**

## Test kits (CE-IVD)

- **GXCARBAR-CE-10**

## Sample collection kits

- **Cepheid catalog number 900-0370**

# Xpert Carba-R Kit

|                                    | Xpert Carba-R Assay   |
|------------------------------------|-----------------------|
| <b>Catalogue Number</b>            | <b>GXCARBAR-CE-10</b> |
| <b>Tests per kit</b>               | 10                    |
| <b>Contents per test cartridge</b> | Dry & Liquid Reagents |
| <b>Transfer pipettes</b>           | 10                    |
| <b>Sample Reagent vials</b>        | 10                    |
| <b>Storage</b>                     | 2-28°C                |





# Xpert Carba-R Sample Collection, Transport, and Storage



Cepheid Sample Collection Device (P/N 900-0370)

SCORE MARK

- **Sample Type:**
  - **Rectal Swab Collection**
- **Sample Collection:**
  - **Collect a paired rectal swab by carefully inserting both swab tips approximately 1 cm beyond the anal sphincter and rotate gently.**
- **Sample Transport and Storage:**
  - **Directly after collection, the swab sample can be held for up to 6 hours at 15– 28°C. Thereafter, the swab can be stored at 2 – 28°C for 7 days.**

# Sample Collection

## Rectal Specimen Collection Protocol


for use with Xpert® assays:

- Xpert Carba-R

GeneXpert®



A Copan Venturi Transystem double swab must be used to collect the specimen.  
Cepheid Collection Device #900-0370



Carefully insert both swab tips approximately 1 cm beyond the anal sphincter and rotate gently.



Place the swabs back in the tube.



Specimens that can be tested within 6 hours can be stored at 15-28 °C. Specimens stored at 2-28 °C are stable for up to 7 days.

defining *on-demand* molecular diagnostics

© 2014 Cepheid

 **Cepheid.**  
Bring answers to life.

301-3343 Rev. A, May 2014

# Good Laboratory Practice

## PCR laboratory setup

- Cartridge/reagent preparation → Sample addition → Detection

## Specimen and reagent storage

- Store specimens separately from reagents to prevent reagent contamination.

## Equipment

- Use filtered pipette tips, when needed, for QC dilutions.
- Follow the manufacturer's recommendation for calibration and maintenance of the lab equipment.

# Good Laboratory Practice, cont'd

## Housekeeping

- Clean work surfaces with a final concentration of 1:10 dilution of household bleach and then a 70% ethanol or 70% isopropanol solution. Wipe work surfaces dry.
- If contamination occurs, thoroughly clean the contaminated area with 1:10 dilution of household bleach, DNA AWAY, or 3% (w/v) hydrogen peroxide and rinse thoroughly with water. Wipe work surfaces dry.

## Personnel

- Wear clean lab coats and gloves.
- Change gloves between processing samples.

## Lab bench area

- Clean the lab bench area routinely.
- Keep the back of the instrument dust free.

# Xpert Carba-R Kit Storage and Handling



- Store test kits at 2-28°C. Do not use expired cartridges.
- Each single-use cartridge is used to process one test. Do not reuse processed cartridges.
- Do not open a cartridge until ready to add the sample eluted from the swab to the cartridge..
  - Cartridge should be placed onto the instrument within 30 minutes of adding the sample into the cartridge.
- Avoid cross contamination during sample handling steps.
  - Change gloves between processing each sample.
  - Change gloves before leaving work area and upon entry into work area.
- Do not use a cartridge that has been dropped after removing it from the packaging.
- Do not shake the cartridge. Shaking or dropping the cartridge after opening the cartridge lid may yield invalid results.
- Do not use a cartridge that has a damaged reaction tube.

# Xpert Carba-R Cartridge Preparation

## Xpert® Cartridge Preparation Carba-R

Refer to the package insert for detailed instructions, precautions, and warnings.

For a copy of the SDS, visit [www.cepheid.com](http://www.cepheid.com) or [www.cepheidinternational.com](http://www.cepheidinternational.com)



Cepheid Technical Support

US office  
(888) 838-3222  
[techsupport@cepheid.com](mailto:techsupport@cepheid.com)

European office  
+33 563 82 53 19  
[support@cepheideurope.com](mailto:support@cepheideurope.com)



- 1 Obtain one Xpert cartridge and one Sample Reagent vial for each sample.
- 2 Insert the swab into the Sample Reagent vial.
- 3 Break the swab at the score mark near the opening of the vial.
- 4 Recap the Sample Reagent vial and vortex for 10 seconds.
- 5 Open the Xpert cartridge lid.
- 6 Aspirate the Sample Reagent up to the line on the supplied pipette.
- 7 Empty the pipette into the sample chamber.
- 8 Close the Xpert cartridge lid.
- 9 Start the test within the timeframe specified in the package insert.




*Note: Do not hold the swab below the score mark. Use gauze or its equivalent to minimize the risk of contamination.*

© 2014 Cepheid

301-3342 Rev. A April 2014

# Running the Test

- Place a single swab into the Sample Reagent vial; break off at score-mark using sterile gauze to minimize risks of contamination.
  - Place the unused swab into the transport tube and store at 2-28°C.
- Vortex the Sample Reagent vial for 10 seconds at high speed.
- Open the Sample Reagent vial.
- Fill the transfer pipette to the mark on the pipette.
  - Avoid air bubbles.
  - The remaining sample in the sample reagent vial can be retained at 2-28°C for up to 4 days in case a retest is required.



# Running the Test- new cartridge image





# Quality Control

Refer to the Package Insert for further details



 **Cepheid**<sup>®</sup>  
*A better way.*

# Cepheid Assay Control Strategy

- **Each Xpert cartridge is a self-contained test device.**
  - **Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge.**
    - Instrument system control: Check status
    - Reagent control: Probe Check
    - Sample processing control: SPC and/or SAC
    - Amplification control: SPC and/or SAC and/or IC

# Instrument System Control – Check Status

- **The Instrument System control checks the optics, temperature of the module, and mechanical integrity of each cartridge.**
  - **If the system controls fail, an ERROR test result will be reported.**

# Reagent Control - Probe Check Control

- After sample preparation, bead reconstitution, and tube filling (prior to thermal cycling), multiple fluorescent readings are taken at different temperatures.
- The readings are compared to default settings established by Cepheid.
- The Probe Check controls for:
  - Missing Target Specific Reagent (TSR) and/or Enzyme Reagent beads, which contain all primers, probes, and internal control template
  - Incomplete reagent reconstitution
  - Incomplete reaction tube fill
  - Probe degradation
- If the Probe Check fails, an ERROR test result will be reported.

# Sample Processing Control - SPC

- The Sample Processing Control (SPC) assesses the effectiveness of the sample preparation steps, including reaction tube filling.
- SPC is *Bacillus globigii* spores.
- The SPC controls for:
  - Missing primer/probe or enzyme beads
  - Incomplete reagent reconstitution
  - Incomplete reaction tube fill
  - Enzyme degradation
  - Sample lysis, nucleic acid extraction, and integrity of nucleic acid
  - Sample inhibition
- The SPC can be negative or positive in an analyte-positive sample.
- If the SPC fails in an analyte-negative sample, an INVALID test result will be reported.

# Commercially Available External Controls

| Organism Name                 | Beta-lactamases Present | Source        |
|-------------------------------|-------------------------|---------------|
| <i>Klebsiella pneumoniae</i>  | KPC-3                   | NCTC 13438    |
| <i>Klebsiella pneumoniae</i>  | KPC                     | ATCC BAA-1705 |
| <i>Escherichia coli</i>       | IMP                     | NCTC 13476    |
| <i>Pseudomonas aeruginosa</i> | VIM-10                  | NCTC 13437    |
| <i>Klebsiella pneumoniae</i>  | VIM-1                   | NCTC 13439    |
| <i>Klebsiella pneumoniae</i>  | VIM-1                   | NCTC 13440    |
| <i>Klebsiella pneumoniae</i>  | NDM-1                   | NCTC 13443    |
| <i>Klebsiella pneumoniae</i>  | NDM-1                   | ATCC BAA-2146 |
| <i>Klebsiella pneumoniae</i>  | OXA-48                  | NCTC 13442    |

- **Other options:**
  - **Known patient positive and negative samples**

# Results Analysis

Refer to the Package Insert for further details




 **Cepheid**<sup>®</sup>  
*A better way.*

# Xpert Carba-R Results:

## IMP1 Pos/ VIM Pos/ NDM Pos/ KPC Pos/ OXA48 Pos

- **IMP1, VIM, NDM, KPC and OXA48 target DNA sequences are detected.**
- **SPC: Not applicable (N/A).** The SPC is ignored because IMP1, VIM, NDM, KPC and OXA48 target DNA amplifications can compete with this control.
- **Probe Check: PASS.** All probe check results pass.






# Xpert Carba-R Results:

## IMP1 Pos/ VIM Neg/ NDM Neg/ KPC Neg/ OXA48 Neg


- **IMP1 target DNA sequence is detected; VIM, NDM, KPC and OXA48 target DNA sequences are not detected.**
- **SPC: Not applicable (NA). SPC is ignored because IMP1 target DNA amplification can compete with this control.**
- **Probe Check: PASS. All probe check results pass.**



# Xpert Carba-R Results:


## IMP1 Pos/ VIM Neg/ NDM Pos/ KPC Neg/ OXA48 Neg

- **IMP1 and NDM target DNA sequences are detected; VIM, KPC and OXA48 target DNA sequences are not detected.**
- **SPC: Not applicable (NA). SPC is ignored because IMP1 and NDM target DNA amplification can compete with this control.**
- **Probe Check: PASS. All probe check results pass.**




# Xpert Carba-R Results: IMP1 Neg/ VIM Neg/ NDM Neg/ KPC Neg/ OXA48 Neg

- **IMP1, VIM, NDM, KPC and OXA48 target DNA sequences are not detected.**
- **SPC: PASS; PCR amplification of the SPC DNA sequence gives a Ct value within the valid range and a fluorescence endpoint above the minimum setting.**
- **Probe Check: PASS. All probe check results pass.**



# Xpert Carba-R Results: INVALID

- **Presence or absence of IMP1, VIM, NDM, KPC and OXA48 target DNA sequences cannot be determined.**
- **SPC: FAIL; No PCR amplification of the SPC DNA sequence or the SPC Ct is not within valid range and the fluorescence endpoint is below minimum setting.**
- **PCC: PASS; all probe check results pass.**



# Xpert Carba-R Result: ERROR

- Presence or absence of IMP1, VIM, NDM, KPC and OXA48 target DNA sequences cannot be determined.
  - SPC: NO RESULT
  - PCC: FAIL\*; all or one of the probe check results fail. The PCC probably failed because the reaction tube was filled improperly or a probe integrity problem was detected.
- \* If the probe check passed, the error is caused by a system component failure.

The screenshot displays the Xpert Carba-R software interface. At the top, there are tabs for 'Test Result', 'Analyte Result', 'Detail', 'Errors', 'History', and 'Support'. Below the tabs, the 'Assay Name' is 'Xpert CARBA-R' and the 'Version' is 'Version 1'. The 'Test Result' field shows 'ERROR' in yellow text. Below this, there is a section labeled 'For In Vitro Diagnostics Use Only'. A 'Troubleshoot' button is visible below the main interface. A table with the following data is shown:

| # | Description          | Detail  | Time                 |
|---|----------------------|---|----------------------|
| 1 | Operation terminated | Error 2008: Syringe pressure reading of 100.0 PSI exceeds the protocol limit of 100.0 PSI | 05/08/13<br>15:51:28 |

## Xpert Carba-R Result: NO RESULT

- Presence or absence of Carba-R target DNA cannot be determined. Use the instructions in Section 13, **Retest Procedure** in the Package Insert to repeat the test. Insufficient data were collected to produce a test result (for example, the operator stopped a test that was in progress).
- SPC: NO RESULT
- PCC: Not applicable

# Reasons to Repeat the Assay

- An **INVALID** result indicates that the SPC control failed. The sample was not properly processed or PCR was inhibited.
- An **ERROR** result indicates that the Probe Check control failed and the assay was aborted possibly due to the reaction tube being filled improperly, a reagent probe integrity problem was detected, or because the maximum pressure limits were exceeded.
- A **NO RESULT** indicates that insufficient data were collected. For example, the operator stopped a test that was in progress.

# Factors That Negatively Affect Results






- **Improper specimen collection**
  - Performance with other collection devices and specimen types has not been assessed.
  - For assays that contain the SAC control, a specimen that does not contain human cells will result in an invalid test result.
- **Improper transport or storage of collected specimen**
  - Storage and transport conditions are specimen specific.
  - Refer to the Package Insert for the appropriate handling instructions.
- **Improper testing procedure**
  - Modification to the testing procedures may alter the performance of the test.
  - Technical error or sample mix-up can impact test results.
  - Careful compliance with the Package Insert is necessary to avoid erroneous results.
- **Interfering substance**
  - False negative test results or invalid results may be observed in the presence of interfering substance.
- **The number of organisms in the specimen is below the detection limit of the test**
- **Refer to Package Insert for non-determinate rate**



# Interfering Substances

- **Of the 23 potentially inhibitory substances tested, Pepto-Bismol 0.25%w/v had a statistically significant inhibitory effect on the detection of IMP-1 in the Xpert Carba-R Assay. No other statistically significant inhibitory effects were observed.**
- **Please refer to the Xpert Carba-R Package Insert for additional data on potentially interfering substances.**

# Carba-R Retest Procedure

|   |   |   |
|---|---|---|
| 1 | Discard used cartridge.   |    |
| 2 | From the test kit, remove new: Xpert Carba-R cartridge, Sample Reagent vial, and transfer pipette.                                  |    |
| 3 | Transfer the remaining liquid from the original Sample Reagent vial (stored $\leq 4$ days at 2-8°C) to the new Sample Reagent vial. |    |
| 4 | Repeat the test with a new cartridge.   |  |
| 5 | Follow the Package Insert on how to run a test.   |  |

## Discussion and Q&A



 **Cepheid<sup>®</sup>**  
*A better way.*

# Technical Support

- **Cepheid provides technical support in the field, on the phone, by fax, and by email.**
- **North America:**
  - **Phone: +1 (888) 838-3222**
  - **Fax: +1 (408) 716-2550**
  - **Email: [techsupport@cepheid.com](mailto:techsupport@cepheid.com)**
- **Europe, Africa, Middle-East, Indian Peninsula:**
  - **Phone:**
    - France: +33 563 825 319**
    - Germany : +49 69 710 480 480**
    - UK: +44 3303 332 533**
    - South Africa: +27 11 234 9636**
    - Australia, New Zealand : +61 1800 107 884**
    - Other European, Middle East and African countries: +971 4 253 3218**
  - **Fax: +33 563 825 301**
  - **Email: [support@cepheideurope.com](mailto:support@cepheideurope.com)**

# Discussion and Q&A

