

# Assay Technical Training

## Xpert<sup>®</sup> HIV-1 Qual

For CE-IVD Use Only

*Cepheid Training Center*



# Training Agenda

- **Xpert HIV-1 Qual**
  - Clinical utility
  - Reagents
  - Sample collection
  - Kit storage and handling
  - Preparing the cartridge
  - Quality Controls
  - Results analysis
- **Discussion**



# Training Objectives

- **At the end of the training, user will be able to:**
  - Properly store and handle the Xpert<sup>®</sup> HIV-1 Qual cartridge kit and sample collection kits
  - Follow proper laboratory safety precautions
  - Collect appropriate specimen types and transport specimen
  - Perform the cartridge set up and run the assay
  - Understand assay control strategy
  - Report the various software generated results

# The Cepheid Solution



- Detects one target
  - HIV-1
- On-board controls for each individual sample
  - Probe Check Control (PCC)
  - Specimen Processing Control (SPC)
  - Sample Volume Adequacy (SVA)
- Results in approximately 90 minutes
- Closed cartridge system minimizes risk of contamination
- On-demand results
- Random access

# Intended Use

- The Xpert HIV-1 Qual Assay, is a **qualitative *in vitro*** diagnostic test designed to detect Human Immunodeficiency Virus Type 1 (HIV-1) total nucleic acids using **human whole blood (WB) and Dried Blood Spots (DBS)** from individuals suspected of HIV-1 infection.
- The Xpert HIV-1 Qual Assay is intended to aid in the diagnosis of HIV-1 infection in conjunction with clinical presentation & other laboratory markers.
- This assay is not intended to be used as a donor screening test for HIV-1 infection.
- This assay uses **reverse transcription polymerase chain reaction (RT-PCR)** technology.

# Targets and Probes


- **Target**

- HIV-1

- **Probes**

- 1 probe binds to the SPC

- 2 probes bind to **3' end** of **5'LTR region** of HIV-1 RNA (do not bind to HIV-2)



- **Depending on the HIV group, one of the probes will bind.**

# Assay Requirements

## GeneXpert Systems

- GeneXpert Dx Software **v4.7b** or higher
- Infinity Dx Software **v6.4b** or higher

## Test Kits (CE-IVD)

- GXHIV-QA-CE-10

## Sample Collection Whole Blood Protocol

- EDTA tubes

## Sample Collection Dried Blood Spots Protocol

DBS Collection Kit (*Whatman 903 Paper card or equivalent*)  
Sterile Scissors and/or forceps  
Eppendorf Thermomixer® C (*Eppendorf reference 5382 000.015*)  
Eppendorf Smartblock™ (vial rack) (*Eppendorf reference 5309 000.007*)

## Other materials

- Personal Protective Equipment (PPE)
- 1:10 Bleach and 70% ethanol or denatured ethanol

## Optional

- Uninterruptible Power Supply /Surge Protector
- Printer

# Good Laboratory Practice

## Personnel Protective Equipment (PPE)

- Wear clean lab coats and gloves
- Change gloves between processing samples

## Lab Bench area

- Clean work surfaces routinely with:
  - ✓ 1:10 dilution of household bleach
  - ✓ 70% Ethanol Solution
- After cleaning, ensure work surfaces are dry

## Specimens, Samples, and Kits Storage

- Store specimens and sample away from kit to prevent contamination

## Equipment(s)

- Use filtered tips when recommended
- Follow the manufacturer's requirements for calibration and maintenance of equipment(s)



# Kit Handling



# Xpert HIV-1 Qual Kit Contents

<b>Catalog Number</b>	<b>GXHIV-QA-CE-10</b>
<b>Tests Per Kit</b>	10
<b>Other reagent</b>	Sample Reagent
<b>Kit CD</b>	Assay Definition File (ADF)
	Assay Import Instructions
	Package Insert (PDF)
<b>Transfer Pipettes</b>	10 (1mL) and 10 (100uL)
<b>Storage</b>	2-28 °C



*Cartridges contain chemically hazardous substances-please see Package Insert and Safety Data Sheet for more detailed information.*

# Xpert HIV-1 Kit Storage and Handling

- Store the Xpert HIV-1 Qual Assay cartridges and reagents at **2–28°C**
- Follow your institution's safety procedures for working with chemicals and handling biological samples
- Do not use Collection Reagent tubes that have not been validated by Cepheid
- Open the Assay cartridge lid only when adding the Sample, close the lid and proceed with processing
- With GeneXpert System, start the test within **30 minutes** after adding the sample to the cartridge
- With the Infinity System, place the cartridge on the conveyor within **30 minutes** after adding the sample.



# Warnings and Precautions

- Do not shake the cartridge
- Do not use a cartridge that... :
  - appears wet, has leaked or if the lid seal appears to have been broken
  - appears damaged
  - has been dropped after removing it from packaging
  - has been dropped or shaken after adding the sample to it
  - has a damaged reaction tube
  - has been used: each cartridge is single-use to process one test
  - is expired
- Do not reuse spent disposable pipettes



# Xpert HIV-1 Qual Limitations

- The performance of Xpert HIV-1 Qual has not been demonstrated from specimens other than the ones validated, i.e. Whole Blood and Dried Blood Spots
- The performance of the Xpert HIV-1 Qual test has not been evaluated with specimens processed by methods other than those described in the package insert.
- Failure to follow assay procedures may lead to false negative results.
- Inhibitors present in the samples may lead to invalid results.

*For detailed information, refer to the current Package Insert*

# Specimen Collection, Storage and Transport



# Whole Blood Collection

- Collect whole blood specimens in EDTA tubes as per manufacturer's instructions
- After blood collection, mix the specimen by gently inverting the tube **10 times**
- A minimum of **100 $\mu$ L** is required for testing



EDTA tube



Heparin tube



# Dried Blood Spot (DBS) Collection

- Collect dried blood spot (DBS) specimen
- DBS should be prepared using Whatman 903, Munktell filter paper cards or equivalent from blood obtained from:
  - a heel-stick
  - finger-stick
  - toe-stick
  - an EDTA-tube
- Ensure that the entire circle is covered with blood (approximately 60–70  $\mu\text{L}$ ). A minimum of two circles should be made from each sample to allow for retesting.
- Air-dry the card at room temperature for a minimum of 4 hours. Package each card in an individual resealable bag with a desiccant sachet in each bag.







<https://www.spotonsciences.com/knowledge-center/dbs-technology/>





For further details about DBS preparation and packaging please refer to related chapter in “*WHO DBS drug resistance testing Manual*”



# Specimen Collection, Transport and Storage

	Prior to testing	<b>Temperature (°C)</b>	<b>Storage Time</b>
	EDTA anticoagulated whole blood		8 hours
			24 hours
			72 hours

Sent to the testing laboratories in individual re-sealable bags containing desiccant

	Prior to testing	<b>Temperature (°C)</b>	<b>Storage Time</b>
	Dried Blood Spots (DBS)	 	12 weeks
			8 weeks

# Cartridge Preparation




# Cartridge Preparation - Whole Blood

1




Label the Sample Reagent (SR) bottle with Sample ID

2



Label the side of the cartridge with the same ID as the SR bottle

3




Open the cartridge lid



4th Mark 1 ml  
3rd Mark 750 µl

*\*HIV-1 Qual Assay  
1 mL Transfer pipette*

4



Use the transfer pipette\* provided (or an automatic pipette) to transfer 750 µL of the sample reagent into the sample chamber of the cartridge

5 x7



Mix the whole blood sample by inverting the vial at least seven times



Upper Bulb

Fill to Here  
100 µl


*\*\*HIV-1 Qual Assay 100 µL  
Transfer Micropipette*

6



Immediately transfer 100 µL of whole blood using the Micropipette\*\* provided (or an automatic pipette) into the same sample chamber

7



Close the lid firmly

Start the test within the time frame specified in the package insert.


# Cartridge Preparation – Dried Blood Spots




Set the Thermomixer C to 56°C




Label the Sample Reagent (SR) bottle with Sample ID




Using sterilized scissors, cut out one DBS from the filter paper card following the dashed line. If the dashed line is perforated, use forceps to detach the DBS



Place one DBS into the SR vial and make sure it is submerged in the buffer



Incubate this vial in the ThermoMixer C for  
•15 minutes at  
•56°C at 500 rpm




Label the side of the cartridge with the same ID as the SR bottle



Open the cartridge lid



*Adapted from 301-5060*



Use the provided transfer pipette\* (or an automatic pipette) to transfer all the liquid from the vial into the sample chamber of the cartridge



Close the lid firmly

Start the test within the time frame specified in the package insert.

# Run a Test

## 1 Create Test

GeneXpert




Start the test within **30** minutes after adding the sample to the cartridge

GeneXpert  
Infinity





Place the cartridge on the conveyor within **30** minutes of adding the sample.

## 2 Scan barcode messages: Cartridge/ Patient and/or Sample ID



*By default, do not click on  
**Manual Entry** or **Cancel***

## 3 Scan the cartridge



*"For complete details on how to run a test, refer to the Package Insert and the GeneXpert Dx or Xpertise Dx Operator Manuals.*

# Create a Test on GeneXpert Dx Software

4 Complete the fields as required

5 The Assay Protocol is selected automatically

6 The module is selected automatically  
**DO NOT CHANGE IT!!!**

7 Click on Start Test

8 A green light will flash on the module  
Load the cartridge into module and close the door

The screenshot shows the 'Create Test' software interface. The fields are as follows:

- Patient ID: [Empty text box]
- Sample ID: [Empty text box]
- Patient ID 2: [Empty text box]
- Last Name: [Empty text box]
- Name: [Empty text box]
- Select Assay: **Xpert HIV-1 Qual** (highlighted with an orange box)
- Select Module: A3 (dropdown menu)
- Reagent Lot ID\*: 16119
- Expiration Date\*: 2016/1/17
- Test Type: Specimen (dropdown menu)
- Sample Type: Other (dropdown menu)
- Notes: [Empty text area]
- Start Test: [Button, highlighted with an orange box and a mouse cursor]
- Scan Cartridge Barcode: [Button]



# Create a Test on Xpertise Dx Software

4 Complete the fields as required

**Order Test - Test Information**


<b>Patient ID</b> patientid	
<b>Sample ID</b> sampleid	
<b>Last Name</b> patient	<b>First Name</b> id
<b>Assay*</b> Xpert HIV-1 Qual	
<b>Reagent Lot ID*</b> 12102	<b>Cartridge S/N*</b> 282769448
<b>Expiration Date*</b> 2018/11/04	<b>Priority</b> Normal
<b>Test Type</b> Specimen	
<b>Sample Type</b> Other	<b>Other Sample Type</b>
<b>Notes</b>	

5 The Assay Protocol is selected automatically

6 Click on SUBMIT



7 Place the cartridge onto the conveyor belt



# Automated Xpert Protocol





# Quality Controls



- **System Control – Check Status**

- System control checks the optics, temperature of the module and mechanical integrity of each cartridge.
- If the system controls fail, an ERROR test result will be reported.

- **Assay Quality Controls**

- Each Xpert cartridge is a self-contained test device
- Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge
  - Sample Volume Adequacy (SVA)
  - Specimen Processing Control (SPC)
  - Probe Check Controls (PCC)

# Internal Quality Controls

- **Sample Volume Adequacy (SVA)**

- Verifies that the correct sample volume is added to the cartridge

- **Probe Check Controls (PCC)**

- Before the PCR step, fluorescence signal is measured on all probes and compared with default factory settings to monitor
  - bead rehydration
  - reaction tube filling
  - probe integrity
  - dye stability

- **Sample Processing Controls (SPC)**

- Armored RNA
  - Verifies adequate sample processing
  - Verifies lysis and presence of the organism and detects PCR inhibition
  - Should be positive in a negative sample
  - Can be positive or negative in a positive sample

# Commercially Available External Controls

## AcroMetrix® <http://www.lifetechnologies.com/acrometrix>

Part Number	Description	Configuration	Storage
964003	HIV-1 High control	1.2 mL x 5 vials	≤ -20°C
964002	HIV-1 Mid control	1.2 mL x 5 vials	≤ -70°C
964001	HIV-1 Low control	1.2 mL x 5 vials	≤ -70°C


1. *Take 1 vial of the control material*
2. *Thaw it at ambient temperature and mix well*
3. *Immediately after thawing place the vial in the ice*
4. *Transfer the entire amount (1.2mL) using the transfer pipette of the Xpert HIV-1 kit and add it to the Xpert HIV-1 sample chamber of the cartridge*
5. *Close the lid and launch the test on GeneXpert*

- *Many other vendors for quality control material are also available than the one outlined above.*
- *External controls should be used in accordance with local, state accrediting organizations, as applicable*

# Result Interpretation



# Result Interpretation Algorithm




# HIV-1 DETECTED

HIV-1 DETECTED

Test Result	Analyte Result	Detail	Errors	History	Support
Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result	
SPC	36.1	476	NA	PASS	
HIV-1	30.8	493	POS	PASS	
QC-1	0.0	0	NEG	PASS	
QC-2	0.0	0	NEG	PASS	

The target HIV-1 is detected


- SPC: NA (Not Applicable) SPC is ignored when HIV-1 target amplified
- Probe Check: PASS



# HIV-1 NOT DETECTED

HIV-1 NOT DETECTED

Test Result	Analyte Result	Detail	Errors	History	Support
Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result	
SPC	34.0	401	PASS	PASS	
HIV-1	0.0	-3	NEG	PASS	
QC-1	0.0	0	NEG	PASS	
QC-2	0.0	0	NEG	PASS	



The target HIV-1 is NOT detected

– SPC: PASS

SPC has a Ct value within the valid range

– Probe Check: PASS



# Troubleshooting




# Factors That Negatively Affect Results

- Improper specimen collection
  - Performance with other collection devices and specimen types has not been assessed
- Improper transport or storage of collected specimen
  - Refer to the Package Insert for the appropriate handling instructions
- Improper testing procedure
  - Modification to the testing procedures, technical error and sample mix-up may impact the test results
  - Careful compliance with the package insert is necessary to avoid erroneous results
- Interfering substance
  - False negative results or invalid results may be observed in the presence of interfering substance

# INVALID Result

INVALID

Test Result	Analyte Result	Detail	Errors	History	Support
Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result	
SPC	39.9	80	FAIL	PASS	
HIV-1	0.0	-3	INVALID	PASS	
QC-1	0.0	0	NEG	PASS	
QC-2	0.0	0	NEG	PASS	



Presence or absence of the HIV-1 target can not be determined

- SPC: FAIL

- SPC Ct value is not within the valid range

- Probe Check: PASS

## - Possible Causes

- Improper sample collection (using heparin tube for e.g)
- Incorrect sample preparation
- Improper storage of the cartridges
- Inefficient sample processing in cartridge
- Presence of inhibitors in the sample

## - Solution

- Repeat the test with a new cartridge and new sample

# Assay Interference

- Potentially Interfering Substances
  - A total of 5 endogenous substances were evaluated
  - Elevated levels of those endogenous substances were shown **not to impact** the assay specificity or interfere with the HIV-1 detection

<b>Substance</b>	<b>Tested Concentration</b>
Albumin (BSA)	90 mg/mL
Bilirubin	0.2 mg/mL
Hemoglobin	5 mg/mL
Human DNA	4 µg/mL
Triglycerides	30 mg/mL

# ERROR Result (2096/2097)

ERROR

Test Result	Analyte Result	Detail	Errors	History	Support
Troubleshoot					
#	Description	Detail			
1	Operation terminated	Error 2097: Assay-Specific Termination Error #2: 46, 29, 1, 0			

The Sample Volume Adequacy (SVA) passes if it meets the validated acceptance criteria.

Error Code	Cause	Solution
<b>2096</b>	No sample added	<ul style="list-style-type: none"><li>– Ensure the Sample is added to cartridge</li><li>– Ensure cartridge is loaded within 30 min. after adding sample</li></ul>
<b>2097</b>	Not enough sample added	<ul style="list-style-type: none"><li>– Ensure the minimum sample volume is added to the cartridge</li><li>– Ensure cartridge is loaded within 30 minutes after adding sample</li></ul>

# NO RESULT

NO RESULT

The screenshot displays a software interface with a top navigation bar containing tabs: Test Result, Analyte Result, Detail, Melt Peaks, Errors, History, Messages, and Support. Below the navigation bar, the 'Assay Name' is 'Xpert HIV-1' and the 'Version' is '1'. The 'Test Result' field is highlighted with a blue border and contains the text 'NO RESULT'. Below this, there is a large empty rectangular area. At the bottom of the interface, the text '<No Data Available>' is displayed.

- The presence or absence of HIV-1 cannot be determined.
  - TARGET 1: NO RESULT
  - CONTROL: NO RESULT
- Probe Check: NA (not applicable)

## Possible Causes

A NO RESULT indicates that insufficient data were collected.

- Test was stopped with stop test button
- Electrical failure

## Solution

- Secure the power
- Repeat the test with a new cartridge

# Re-test Procedure

1



Discard used cartridge

*Follow your institution's safety guidelines for disposal of cartridges*

2



Obtain the residual sample, mix according to Package Insert

*If the leftover sample volume is insufficient, or the retest continues to return an INVALID, ERROR, or NO RESULT, collect a new sample*

3



Obtain a new cartridge


Label appropriately as retest on the new cartridge

Process the sample per the package insert

4



Run the test on the GeneXpert System



# Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
  - Product name
  - Lot number
  - Serial number of the System
  - Error messages (if any)
  - Software version and, if applicable, Computer Service Tag number
- Log you complaint on the following link <http://www.cepheid.com/us/support>

Region	Telephone	Technical Support Email
US	+ 1 888 838 3222	<a href="mailto:techsupport@cepheid.com">techsupport@cepheid.com</a>
Australia and New Zealand	+ 1800 130 821 + 0800 001 028	<a href="mailto:techsupportanz@cepheid.com">techsupportanz@cepheid.com</a>
Brazil and Latin America	+ 55 11 3524 8373	<a href="mailto:latamsupport@cepheid.com">latamsupport@cepheid.com</a>
China	+ 86 400 821 0728	<a href="mailto:techsupportchina@cepheid.com">techsupportchina@cepheid.com</a>
Japan	0120 95 4886	<a href="mailto:support@japan.cepheid.com">support@japan.cepheid.com</a>
France	+ 33 563 825 319	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Germany	+ 49 69 710 480 480	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
India, Bangladesh, Bhutan, Nepal and Sri Lanka	+ 91 11 48353010	<a href="mailto:techsupportindia@cepheid.com">techsupportindia@cepheid.com</a>
Italy	+ 39 800 902 567	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
South Africa	+ 27 861 22 76 35	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
United Kingdom	+ 44 3303 332 533	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Other European, Middle East, and African countries	+ 33 563 825 319 + 971 4 253 3218	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Other countries not listed above	+1 408 400 8495	<a href="mailto:techsupport@cepheid.com">techsupport@cepheid.com</a>





Thank You.



[www.Cepheid.com](http://www.Cepheid.com)

