

# Assay Technical Training

## Xpert<sup>®</sup> Xpress Flu/RSV

For CE-IVD Use Only

*Cepheid Training Center*





# Training Agenda

- **Xpert Xpress Flu/RSV Training**
  - Reagents
  - Sample collection
  - Kit storage and handling
  - Preparing the cartridge
  - Quality Controls
  - Results analysis
- **Discussion**





# Training Objectives

- **At the end of the training, users will be able to:**
  - Store and handle the Xpert<sup>®</sup> Xpress Flu/RSV kit
  - Follow proper laboratory safety precautions
  - Collect and transport appropriate specimen(s)
  - Prepare a cartridge and run the assay
  - Report and understand various software-generated results
  - Understand the assay control strategy



# The Cepheid Solution



- Detection and differentiation of influenza A, influenza B, and respiratory syncytial virus (RSV)
- On-board internal controls for each sample
  - Probe Check Control (PCC)
  - Sample Processing Control (SPC)
- Closed cartridge system minimizes risk of contamination
- On-demand results
- Random access

# Intended Use

- The Cepheid Xpert® Xpress Flu/RSV test, performed on the GeneXpert® Instrument Systems, is an automated, multiplex real-time, reverse transcriptase polymerase chain reaction (RT-PCR) assay intended for the *in vitro* qualitative detection and differentiation of influenza A, influenza B, and respiratory syncytial virus (RSV). The Xpert Xpress Flu/RSV test uses nasopharyngeal (NP) swab or nasal swab (NS) specimens collected from patients with signs and symptoms of respiratory infection. The Xpert Xpress Flu/RSV test is intended as an aid in the diagnosis of influenza and respiratory syncytial virus infections in conjunction with clinical and epidemiological risk factors.
- Negative results do not preclude influenza virus or RSV infection and should not be used as the sole basis for treatment or other patient management decisions.
- Performance characteristics for influenza A were established during the 2015-2016 influenza season. When other novel influenza A viruses are emerging, performance characteristics may vary.
- If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to state or local health departments for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.



# Xpert Xpress Flu/RSV Requirements

## GeneXpert Systems

- GeneXpert Dx Software v **4.7b** or higher
- Xpertise Software v **6.4b** or higher

## Test Kits (CE-IVD)

- XPRSFLU/RSV-CE-10

## Sample Collection

- Xpert Nasopharyngeal Sample Collection Kit (#SWAB/B-100)
- Xpert Swab Sample collection Kit (#SWAB/F-100)

## Other materials

- Personal Protective Equipment (PPE)
- 1:10 bleach
- 70% ethanol or denatured ethanol

## Optional

- Uninterruptible Power Supply /Surge Protector
- Printer
- Vortex



# Good Laboratory Practice

## Personnel Protective Equipment (PPE)

- Wear clean lab coats and gloves
- Change gloves between processing samples

## Lab Bench Area

- Clean work surfaces routinely with:
  - ✓ 1:10 dilution of household bleach
  - ✓ 70% ethanol solution
- \* Final Active Chlorine concentration should be 0.5% regardless of the household bleach concentration in your country
- After cleaning, ensure that the work surfaces are dry

## Specimens, Samples, and Kits Storage

- Store specimens and sample away from kit to prevent contamination

## Equipment

- Use filtered tips when recommended
- Follow the manufacturer's requirements for calibration and maintenance of equipment

# Kit Handling





# Xpert Flu/RSV Kit Contents

Catalog Number	XPRSFLU/RSV-CE-10
Cartridges Per Kit	10
Transfer Pipettes	1 bag of 12 (300µl volume)
Kit CD	Assay Definition File (ADF)
	Assay Import Instructions
	Package Insert (PDF)
Storage	2-28 °C



*Cartridges contain chemically hazardous substances-please see Package Insert and Safety Data Sheet for more detailed information.*



# Xpert Xpress Flu/RSV Kit Storage and Handling

- Store the Xpert Assay cartridges and reagents at 2–28°C
- Follow your institution's safety procedures for working with chemicals and handling biological samples
- Do not use sample collection devices that have not been validated by Cepheid
- Open the cartridge lid only when adding the sample, close the lid, and proceed with processing





# Warnings and Precautions

- Do not shake the cartridge
- Do not use a cartridge that... :
  - appears wet, has leaked, or if the lid seal appears to have been broken
  - appears damaged
  - has been dropped after removing it from packaging
  - has been dropped or shaken after adding the sample to it
  - has a damaged reaction tube
  - has been used; each cartridge is single-use to process one test
  - is expired
- Do not reuse disposable pipettes

# Warnings and Precautions

- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents and require use of standard precautions.
- Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures.
- If national or regional regulations do not provide clear direction on proper disposal, biological specimens and used cartridges should be disposed per WHO [World Health Organization] medical waste handling and disposal guidelines.



# Specimen Collection, Storage and Transport



# Specimen Collection Device

	Nasopharyngeal Swab and Transport Medium	Nasal Swab and Transport Medium
Catalog #	SWAB/B-100	SWAB/F-100
Intended Use	Designed to collect, preserve and transport respiratory virus specimens	
Kit Contents	1 tube swab transport medium (3 mL, red cap)	
	1 flocked swab with flexible shaft	1 flocked swab with rigid shaft




# Specimen Collection- Nasal Swab


## Nasal Swab Specimen Collection

For use with Xpert® Swab Sample Collection Kit - Catalog # SWAB/F-100


- 1 Open the package that contains the swab and transport medium tube. Set the tube aside before collecting the specimen.




5 Repeat Step 4 on the other nostril with the same swab.





6 To avoid specimen contamination, do not touch the swab tip to anything after collecting the specimen.


- 2 Open the swab wrapper and remove the swab, taking care not to touch the tip of the swab to any surface.





7 Remove the cap from the tube. Insert the swab into the transport medium.


- 3 Hold the swab in your hand, pinching in the middle of the swab shaft on the scoreline.




8 Break the swab shaft against the side of the tube at the scoreline.


- 4 Rotate swab against the inside of the nostril for 3 seconds while applying pressure with a finger to the outside of the nostril.



8 Avoid splashing contents on the skin. Wash with soap and water if exposed.



\* SWAB/F-100 contains Copan UTM 330C and Copan nylon swab 502CS01

© 2017 Cepheid

In Vitro Diagnostic Use



In Vitro Diagnostic Use



301-9007, Rev. A January 2016


  
Cepheid.  
A better way.


# Specimen Collection- Nasopharyngeal Swab


## Nasopharyngeal Specimen Collection


For use with Xpert® Nasopharyngeal Sample Collection Kit - Catalog # SWAB/B-100


- 1 ..... Open the package that contains the swab and transport medium tube. Set the tube aside before collecting the specimen.



- 2 ..... Open the swab wrapper and remove the swab, taking care not to touch the tip of the swab to any surface.



- 3 ..... Hold the swab in your hand, pinching in the middle of the swab shaft on the scoreline.


- 4 ..... Gently insert the swab into the nostril until you touch the posterior nasopharynx. Rotate swab several times.



- 6 ..... Remove the cap from the tube. Insert the swab into the transport medium.


- 7 ..... Break the swab shaft against the side of the tube at the scoreline. Avoid splashing contents on the skin. Wash with soap and water if exposed.


- 8 ..... Replace the cap on the tube and close tightly. Specimen should be transported at 2-8 °C. Prior to testing, specimen may be stored for 24 hours at 15-30 °C or up to 7 days at 2-8 °C



\* SWAB/B-100 contains Copan UTM 330C and Copan nylon swab 503CS01

© 2017 Cepheid In Vitro Diagnostic Use 




In Vitro Diagnostic Use  

301-6052, Rev. C, January 2018

  
Cepheid.  
A better way.



# Specimen Collection, Transport and Storage

Sample type	Transportation	Storage Conditions
<b>Xpert Viral Transport Medium containing:  Nasopharyngeal Swab Or Nasal Swab</b>		 Up to 24 hours   Up to 7 days

# Cartridge Preparation



# Xpert Xpress Flu/RSV Cartridge Preparation

## Xpert® Xpress Flu and Flu/RSV Cartridge Preparation

- Xpert® Xpress Flu/RSV
- Xpert® Xpress Flu\*

Refer to the package insert for detailed instructions, precautions, and warnings.

For a copy of the SDS, visit [www.cepheid.com](http://www.cepheid.com) or [www.cepheidinternational.com](http://www.cepheidinternational.com)

Cepheid Technical Support  
US office  
(888) 838-3222, Option 2  
[techsupport@cepheid.com](mailto:techsupport@cepheid.com)

European office  
+33 563 825 319  
[support@cepheid.com](mailto:support@cepheid.com)



- 1 Take one Xpert cartridge for each sample.



- 2 Invert the tube 5 times.



- 3 Open the cartridge lid.



- 4 Using a clean 300  $\mu$ L pipette (supplied), transfer 300  $\mu$ L (one draw), of the sample to the opening of the cartridge.



- 5 Close the cartridge lid.



- 6 Start the test within the timeframe specified in the package insert.

© 2017 Cepheid


\*FDA Cleared Assay Only

301-7307, Rev. B February 2017


# Xpress Flu/RSV Cartridge Preparation



Take one Xpert cartridge for each sample.




Invert the tube 5 times.



Open the cartridge lid.



Using a clean 300  $\mu$ L pipette (supplied), transfer 300  $\mu$ L (one draw), of the to the sample to the cartridge.



Close the cartridge lid.

6

Start the test within the timeframe specified in the package insert.



# Run a Test

## 1 Create Test

GeneXpert




Start the test within **30 minutes** after adding the sample to the cartridge

GeneXpert  
Infinity



Place the cartridge on the conveyor within **30 minutes** of adding the sample.

## 2 Scan barcode : Cartridge/ Patient and/or Sample ID



*By default, do not click on  
Manual Entry or Cancel*

## 3 Scan the cartridge



*"For complete details on how to run a test, refer to the Package Insert and the GeneXpert Dx or Xpertise Operator Manuals.*

# Create a Test on GeneXpert Dx Software

4 Complete the fields as required

5 The Assay Protocol is selected automatically

6 The module is selected automatically

7 Click on Start Test

8 A green light will flash on the module  
Load the cartridge into module and close the door

The screenshot shows the 'Create Test' window with the following fields and values:

- Patient ID: [Empty]
- Sample ID: [Empty]
- Patient ID 2: [Empty]
- Last Name: [Empty]
- Name: [Empty]
- Select Assay: Xpert Assay name
- Select Module: A3
- Reagent Lot ID\*: 16119
- Expiration Date\*: 2016/1/17
- Test Type: Specimen
- Sample Type: Other
- Notes: [Empty]

At the bottom, there are two buttons: 'Start Test' and 'Scan Cartridge Barcode'. A mouse cursor is pointing at the 'Start Test' button.



# Combinatorial reporting: Xpert Xpress Flu/RSV

- Choose the desired test from the “Select Assay” drop-down menu

Name	Version
Xpert Xpress Flu-RSV	1
Xpert Xpress_RSV	1
Xpert Xpress_Flu	1
Xpert Xpress Flu-RSV	1

- Only the test result for the assay selected at this step will be collected once the test is started.
  - Example: If the operator selects Xpert Xpress\_RSV, once the assay starts, the option cannot be changed to collect Flu data

# Create a Test on Xpertise Software

4 Complete the fields as required


**Order Test - Test Information**

<b>Patient ID</b> patientid	
<b>Sample ID</b> sampleid	
<b>Last Name</b> patient	<b>First Name</b> id
<b>Assay*</b> Xpert Assay	
<b>Reagent Lot ID*</b> 12102	<b>Cartridge S/N*</b> 282769448
<b>Expiration Date*</b> 2018/11/04	<b>Priority</b> Normal
<b>Test Type</b> Specimen	
<b>Sample Type</b> Other	<b>Other Sample Type</b> 
<b>Notes</b>	

5 The assay protocol is selected automatically

6 Click on SUBMIT

7 Place the cartridge into the conveyor belt



# Automated Xpert Xpress Flu/RSV Test Steps





# Quality Controls





# Assay Control Strategy

CONTROL

## Xpert Assay Quality Controls

- Each Xpert cartridge is a self-contained test device
- Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge
  - Probe Check Controls (PCC)
  - Sample Processing Control (SPC)

*Refer to 301-4868 GeneXpert Quality Control Features for All Cepheid Xpert Assays*



# Internal Quality Controls

## Probe Check Controls (PCC)

- Before the PCR step, fluorescence signal is measured on all probes and compared with default factory settings to monitor
  - bead rehydration
  - reaction tube filling
  - probe integrity
  - dye stability

## Sample Processing Control (SPC)

- Verifies that conditions for adequate amplification process were met
- Detects PCR inhibition
- Should be positive in a negative sample
- Can be positive or negative in a positive sample

# Commercially Available External Controls

Zeptomatrix Part Number	Description	Configuration
NATFLUA/B-6C	FLU A/B Positive	6 X 0.5 mL/box
NATCXVA9-6C	Negative Control	6 X 0.5 mL/box
NATRSV-6C	RSV Positive	6 X 0.5 mL/box
<a href="http://www.zeptomatrix.com">http://www.zeptomatrix.com</a>		

- Storage conditions for external controls: 2-8°C
- Other options: **known positive and negative patient samples**

# Result Interpretation





# Results Summary

Result displayed	Flu A 1	Flu A 2	Flu B	RSV	SPC
Flu A POSITIVE; Flu B NEGATIVE; RSV NEGATIVE	+	+/-	-	-	+/-
	+/-	+			
Flu A POSITIVE; Flu B POSITIVE; RSV NEGATIVE	+	+/-	+	-	+/-
	+/-	+			
Flu A POSITIVE; Flu B NEGATIVE; RSV POSITIVE	+	+/-	-	+	+/-
	+/-	+			
Flu A POSITIVE; Flu B POSITIVE; RSV POSITIVE	+	+/-	+	+	+/-
	+/-	+			
Flu A NEGATIVE; Flu B POSITIVE; RSV NEGATIVE	-	-	+	-	+/-


# Results Summary continued

Result displayed	Flu A 1	Flu A 2	Flu B	RSV	SPC
Flu A NEGATIVE; Flu B NEGATIVE; RSV POSITIVE	-	-	-	+	+/-
Flu A NEGATIVE; Flu B POSITIVE; RSV POSITIVE	-	-	+	+	+/-
Flu A NEGATIVE; Flu B NEGATIVE; RSV NEGATIVE	-	-	-	-	+
INVALID	-	-	-	-	-
ERROR	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT
NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT

# Flu A POSITIVE

Test Result  
Flu A POSITIVE;  
Flu B NEGATIVE;  
RSV NEGATIVE


- Flu A target RNA detected;  
Flu B target RNA not detected;  
RSV target RNA not detected
- Flu A POSITIVE; Flu A target has a valid Ct.
- SPC Not Applicable; The SPC is ignored because the Flu A target amplification may compete with this control.
- PCC PASS; All probe check results pass



# Flu B POSITIVE

Test Result  
Flu A NEGATIVE;  
Flu B POSITIVE;  
RSV NEGATIVE


- Flu B target RNA detected;  
Flu A target RNA not detected;  
RSV target RNA not detected
- Flu B POSITIVE; Flu B target has a valid Ct.
- SPC Not Applicable ;The SPC is ignored because the Flu B target amplification may compete with this control.
- PCC PASS; All probe check results pass



# RSV POSITIVE

Test Result  
Flu A NEGATIVE  
Flu B NEGATIVE  
RSV POSITIVE

- RSV target RNA detected; Flu A and Flu B target RNA not detected
- RSV POSITIVE; RSV target has a valid Ct.
- SPC Not Applicable; The SPC is ignored because the RSV target amplification may compete with this control.
- PCC PASS; All probe check results pass






# Flu A and Flu B POSITIVE

Test Result	Flu A POSITIVE; Flu B POSITIVE; RSV NEGATIVE
-------------	--


- Flu A and Flu B target RNA detected;  
RSV target RNA not detected
- Flu A and Flu B targets have a valid Ct.
- SPC Not Applicable; The SPC is ignored because the Flu A and Flu B target amplification may compete with this control.
- PCC PASS; All probe check results pass



# Flu A and RSV POSITIVE

Test Result  
Flu A POSITIVE;  
Flu B NEGATIVE;  
RSV POSITIVE


- Flu A and RSV target RNA detected;  
Flu B target RNA not detected
- Flu A and RSV targets have a valid Ct.
- SPC Not Applicable; The SPC is ignored because the Flu A and RSV target amplification may compete with this control.
- PCC PASS; All probe check results pass



# Flu B and RSV POSITIVE

Test Result  
Flu A NEGATIVE;  
Flu B POSITIVE;  
RSV POSITIVE


- Flu B and RSV target RNA detected;  
Flu A target RNA not detected
- Flu B and RSV targets have a valid Ct.
- SPC Not Applicable; The SPC is ignored because the Flu B and RSV target amplification may compete with this control.
- PCC PASS; All probe check results pass



# Flu A, Flu B, and RSV POSITIVE

Test Result  
Flu A POSITIVE;  
Flu B POSITIVE;  
RSV POSITIVE


- Flu A, Flu B, and RSV target RNA detected
- Flu A, Flu B, and RSV targets have a valid Ct.
- SPC Not Applicable; The SPC is ignored because the Flu A, Flu B, and RSV target amplification may compete with this control.
- PCC PASS; All probe check results pass



# Flu A, Flu B, and RSV Negative

Test Result  
Flu A NEGATIVE;  
Flu B NEGATIVE;  
RSV NEGATIVE

- Flu A, Flu B, and RSV target RNA not detected
- SPC PASS; SPC has a Ct within the valid range and endpoint above the threshold setting.
- PCC PASS; All probe check results pass





# Xpert<sup>®</sup> Xpress Flu/RSV

## *Early Assay Termination (EAT)*

- EAT for Flu-only and RSV-only ADFs
  - positive results available as early as 20 mins depending on viral titer
- EAT activated when pre-determined threshold for a positive test result is reached before the full 40 PCR cycles completed


### Please note:

- When Flu A or Flu B titers are very high, generating very early Cts with Xpert Xpress Flu, SPC amplification curves may not be seen and not reported.
- When RSV titers are very high, generating very early Cts with Xpert Xpress RSV, SPC amplification curves may not be seen and not reported.

# Flu A POSITIVE and Flu B NEGATIVE

Test Result  
Flu A POSITIVE;  
Flu B NEGATIVE


- Flu A target RNA detected.
- Flu A targets have a valid Ct.
- SPC Not Applicable; The SPC is ignored because the Flu A and Flu B target amplification may compete with this control.
- PCC PASS; All probe check results pass



# RSV POSITIVE

Test Result **RSV POSITIVE**

- RSV target RNA detected
- RSV POSITIVE; RSV target has a valid Ct.
- SPC Not Applicable; The SPC is ignored because the RSV target amplification may compete with this control.
- PCC PASS; All probe check results pass



# Troubleshooting




# Reasons to Repeat the Assay

- An **INVALID** result indicates that the sample was not processed properly, PCR was inhibited, or the sample was inadequate.
- An **ERROR** result indicates that the Probe Check Control failed or maximum pressure limits were exceeded.
- A **NO RESULT** indicates that insufficient data were collected. For example, the operator stopped a test that was in progress, a load error occurred, or the software was closed prematurely.
- Because the incidence of co-infection with two or more viruses (Influenza A, Influenza B, and RSV) is low, it is recommended that specimens undergo repeat testing if nucleic acids from two or more analytes are detected in a single specimen.

# INVALID

Test Result **INVALID**

- SPC does not meet acceptance criteria. Presence or absence of the target RNAs cannot be determined. Repeat the test according to the instructions in the Retest Procedure section of the package insert.
- SPC FAIL; SPC target result is negative. The SPC Ct is not within valid range and fluorescence endpoint is below the minimum setting.
- PCC PASS; all probe check results pass.





# ERROR

Test Result **ERROR**

Presence or absence of Flu A, Flu B, and/or RSV target RNAs cannot be determined.

Repeat the test according to the instructions in the Retest Procedure section of the package insert.

- Flu A NO RESULT
- Flu B NO RESULT
- RSV NO RESULT
  
- SPC NO RESULT
  
- PCC FAIL.\* All or one of the probe check results fail.

\* If the probe check passed, the error is caused by the maximum pressure limit exceeding the acceptable range or by a system component failure.



# NO RESULT

Test Result **NO RESULT**

- The presence or absence of Flu A, Flu B, and RSV target RNAs cannot be determined.
- Repeat the test according to the instructions in the Retest Procedure section of the package insert.
- Flu A NO RESULT
- Flu B NO RESULT
- RSV NO RESULT
- SPC NO RESULT
- PCC Not Applicable

The image shows a screenshot of a diagnostic device's display. At the top, there is a header area with the text "Test Result" followed by a box containing "NO RESULT". Below this, a larger box contains the text "For In Vitro Diagnostics Use Only." The main body of the display is mostly blank, with the text "<No Data Available>" centered at the bottom.

# Re-test Procedure

1

Discard used cartridge

Follow your institution's safety guidelines for disposal of cartridges

2



Obtain the residual sample, mix according to Package Insert

If the leftover sample volume is insufficient, or the retest continues to return an INVALID, ERROR, or NO RESULT, collect a new sample

3



Obtain a new cartridge


Label appropriately as retest on the new cartridge

Process the sample per the Package Insert

4



Run the test on the System





# Factors That Negatively Affect Results

- Improper specimen collection
  - The viral load in the specimen is below the detection limit of the test
  - Performance with other specimen types has not been assessed
- Improper transport or storage of collected specimen
  - Storage and transport conditions are specimen specific
  - Refer to the Package Insert for the appropriate handling instructions
- Improper testing procedure
  - Modification to the testing procedures may alter the performance of the test
  - Careful compliance with the Package Insert is necessary to avoid erroneous results



# Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
  - Product name
  - Lot number
  - Serial number of the System
  - Error messages (if any)
  - Software version and, if applicable, Computer Service Tag number
- Log your complaint online using the following link <http://www.cepheid.com/us/support> : *Create a Support Case*

Region	Telephone	Technical Support Email
US	+ 1 888 838 3222	<a href="mailto:techsupport@cepheid.com">techsupport@cepheid.com</a>
Australia and New Zealand	+ 1800 130 821 (AU) + 0800 001 028 (NZ)	<a href="mailto:techsupportANZ@cepheid.com">techsupportANZ@cepheid.com</a>
Brazil and Latin America	+ 55 11 3524 8373	<a href="mailto:latamsupport@cepheid.com">latamsupport@cepheid.com</a>
China	+ 86 400 821 0728	<a href="mailto:techsupportchina@cepheid.com">techsupportchina@cepheid.com</a>
France	+ 33 563 825 319	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Germany	+ 49 69 710 480 480	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
India, Bangladesh, Bhutan, Nepal, and Sri Lanka	+ 91 1148353010	<a href="mailto:techsupportindia@cepheid.com">techsupportindia@cepheid.com</a>
Italy	+ 39 800 902 567	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
South Africa	+ 27 861 22 76 35	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
United Kingdom	+ 44 3303 332 533	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Belgium and Netherlands	+33 563 825 319	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Other European, Middle East, and African countries	+ 33 563 825 319 + 971 4 253 3218	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Other countries not listed	+ 1 408 400 8495	<a href="mailto:techsupport@cepheid.com">techsupport@cepheid.com</a>



Thank You.



[www.Cepheid.com](http://www.Cepheid.com)

