

Assay Training: Xpert[®] EV

For CE-IVD and US-IVD Use Only



Training Agenda

- Xpert EV Training
 - Reagents
 - Sample collection
 - Kit storage and handling
 - Preparing the cartridge
 - Quality control
 - Results analysis
- Discussion



Xpert EV Training Objectives

At the end of the training, users will be able to:

- Store and handle the Xpert EV cartridge kit.
- Follow proper laboratory safety precautions.
- Identify appropriate specimen types and transport specimens.
- Prepare a cartridge and run the assay.
- Report and understand the various software generated results.
- Understand the assay control strategy.

Xpert® EV



The Cepheid Solution



- Detection of enterovirus (EV) RNA in cerebrospinal fluid (CSF) specimens
- On-board internal controls for each sample
 - Probe Check Control (PCC)
 - Sample Processing Control (SPC)
- Closed cartridge system minimizes risk of contamination
- On-demand results
- Random access

Intended Use

The Cepheid Xpert EV assay is a reverse transcription polymerase chain reaction (RT-PCR) using the GeneXpert® Dx System for the **presumptive qualitative detection of enterovirus (EV) RNA in cerebrospinal fluid (CSF) specimens from individuals with signs and symptoms of meningitis**. This test, in conjunction with other laboratory results and clinical information, may be used as an aid in the laboratory diagnosis of enterovirus infection in patients with a clinical suspicion of meningitis or meningoencephalitis. Assay performance characteristics have not been established for immunocompromised or immunosuppressed patients.

CAUTION: The results obtained with the Xpert EV assay should be used only as an adjunct to clinical observations and other information available to the physician. Positive Xpert EV results do not rule out other causes of meningitis, including bacteria, mycobacteria, other viruses (e.g., herpes family viruses, arboviruses, mumps virus, etc.), and fungi.

System and Reagent Requirements

GeneXpert Systems

- GeneXpert Software v2.1 or higher

Test Kits (US-IVD and CE-IVD)

- GXEV-100N-10

Materials Required but not Provided

- 200- μ L pipette
- Sterile 200- μ L filter barrier pipette tips
- Personal Protective Equipment (PPE)
- 1:10 dilution of bleach
- 70% ethanol or denatured ethanol

Optional

- Uninterruptible Power Supply/ Surge Protector
- Printer

Good Laboratory Practice

Personnel Protective Equipment (PPE)

- Wear clean lab coats and gloves
- Change gloves between processing samples

Lab Bench area

- Clean work surfaces routinely with:
 - ✓ 1:10 dilution of household bleach
 - ✓ 70% Ethanol Solution
- After cleaning, ensure work surfaces are dry

Specimens, Samples, and Kits Storage

- Store specimens and sample away from kit to prevent contamination

Equipment(s)

- Use filtered tips when recommended
- Follow the manufacturer's requirements for calibration and maintenance of equipment(s)

Kit Handling



Xpert EV Kit Contents

Xpert EV

Catalog Number	GXEV-100N-10
Cartridges per Kit	10
Reagent Vials	Binding Reagent
	Wash Reagent
	Elution Reagent
	Lysis Reagent
Kit CD	Assay Definition File (ADF)
	Assay Import Instructions
	Package Insert (PDF)
Storage	2- 28 °C



Xpert EV Kit Storage and Handling

- Store the Xpert EV cartridges and reagents at 2–28°C
- Follow your institution's safety procedures for working with chemicals and handling biological samples
- Open the cartridge lid only when adding the sample, close the lid, and proceed with processing

Warnings and Precautions

- Do not shake the cartridge
- Do not use a cartridge... :
 - if it appears wet, has leaked, or if the lid seal appears to have been broken
 - if it appears damaged
 - that has been dropped after removing it from packaging
 - that has been dropped or shaken after you have added the sample
 - that has a damaged reaction tube
 - that has been used; each cartridge is single-use to process one test
 - is expired
- Do not reuse pipettes



Waste Disposal

- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents and require use of standard precautions.
- Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures.
- If national or regional regulations do not provide clear direction on proper disposal, biological specimens and used cartridges should be disposed per WHO [World Health Organization] medical waste handling and disposal guidelines.

Specimen Storage and Transport



Specimen Handling

Specimen	Storage
Cerebrospinal fluid (CSF) in sterile container	2-8 °C for up to 72 hours following specimen collection
	The specimen may also be frozen at -20 or -80 °C if test will not be performed within 72 hours *Do not freeze and thaw specimen more than twice



Centrifugation of the specimen is not recommended.

Xpert EV Testing Protocol

Xpert EV Cartridge Preparation

Refer to the package insert for detailed instructions, precautions, and warnings.

For a copy of the MSDS, visit www.cephheid.com or www.cephheidinternational.com

Cepheid Technical Support

US office
(888) 838-3222, Option 2
techsupport@cephheid.com

European office
+33 563 82 53 19
support@cephheid.com




- 1 Obtain one cartridge, three ampoules, and one vial for each sample.




- 2 Open the cartridge lid.



- 3 Open, and then add ampoule 1 into 1, ampoule 2 into 2, and ampoule 3 into 3.



- 4 Pipette 140 μ L of the Lysis Reagent into chamber 4S, and then pipette 140 μ L of the sample into chamber 4S.



- 5 Close the cartridge lid.



- 6 Start the test within the timeframe specified in the package insert.


EV Cartridge Preparation




Obtain one cartridge, three ampoules, and one vial for each sample.



Open the cartridge lid.



Open, and then add ampoule 1 into 1, ampoule 2 into 2, and ampoule 3 into 3.



Pipette 140 μ L of the Lysis Reagent into chamber 4S, and then pipette 140 μ L of the sample into chamber 4S.



Close the cartridge lid.

6

Start the test within the timeframe specified in the package insert.

Run a Test

1 Create Test

GeneXpert




Start the test within **30 minutes** after adding the sample to the cartridge

GeneXpert
Infinity





Place the cartridge on the conveyor within **30 minutes** of adding the sample.

2 Scan barcodes: Cartridge/ Patient and/or Sample ID



*By default, do not click on
Manual Entry or Cancel*

3 Scan the cartridge



"For complete details on how to run a test, refer to the Package Insert and the GeneXpert Dx or Xpertise Operator Manuals.

Create a Test on GeneXpert Dx Software

4 Complete the fields as required

5 The Assay Protocol is selected automatically

6 The module is selected automatically

7 Click on Start Test

8 A green light will flash on the module
Load the cartridge into module and close the door

The screenshot shows the 'Create Test' software interface. The fields are as follows:

- Patient ID: [Empty text box]
- Sample ID: [Empty text box]
- Patient ID 2: [Empty text box]
- Last Name: [Empty text box]
- Name: [Empty text box]
- Select Assay: **Xpert Assay name** (highlighted with an orange box)
- Select Module: A3 (dropdown menu)
- Reagent Lot ID*: 16119
- Expiration Date*: 2016/1/17
- Test Type: Specimen (dropdown menu)
- Sample Type: Other (dropdown menu)
- Other S: [Empty text box]
- Notes: [Empty text area]
- Start Test: [Button, highlighted with an orange box and a mouse cursor]
- Scan Cartridge Barco: [Button]



Create a Test on Xpertise Software

4 Complete the fields as required

Order Test - Test Information


Patient ID patientid	
Sample ID sampleid	
Last Name patient	First Name id
Assay* Xpert Assay	
Reagent Lot ID* 12102	Cartridge S/N* 282769448
Expiration Date* 2018/11/04	Priority Normal
Test Type Specimen	
Sample Type Other	Other Sample Type
Notes	



5 The Assay Name Protocol is selected automatically

6 Click on SUBMIT

7 Place the cartridge into the conveyor belt



Automated Xpert Protocol



Quality Controls

*Refer to the Package Insert for
complete details*



Cepheid Assay Control Strategy

CONTROL

- Each Xpert cartridge is a self-contained test device.
- Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge.
 - Probe Check Control (PCC)
 - Sample Processing Control (SPC)

- **Probe Check Controls (PCC)**

- Before the PCR step, fluorescence signal is measured on all probes and compared with default factory settings to monitor
 - bead rehydration
 - reaction tube filling
 - probe integrity
 - dye stability

- **Sample Processing Controls (SPC)-** displayed as CIC

- Verifies adequate sample processing
- Verifies lysis, presence of the organism and detects PCR inhibition
- Should be positive in a negative sample
- Can be positive or negative in a positive sample

Results Analysis

*Refer to the Package Insert for
complete details*



Results Summary

Result displayed	EV	CIC
POSITIVE	+	+/-
NEGATIVE	-	+
INVALID	-	-
ERROR	NO RESULT	NO RESULT
NO RESULT	NO RESULT	NO RESULT


EV Positive

Test Result **POSITIVE**

EV target nucleic acid is detected.

- CIC (SPC/IC)—NA
- Probe Check—PASS
All probe check results pass.

Positive Xpert EV results do not rule out other causes of meningitis, including bacteria, mycobacteria, other viruses, and fungi.




EV Negative

Test Result **NEGATIVE**

EV target nucleic acid is not detected.
EV—NEG

- CIC (SPC/IC) – PASS; SPC has a Ct within the valid range and endpoint above the minimum setting.
- Probe Check – PASS; all probe check results pass.

Negative Xpert EV results do not rule out enterovirus as the cause of meningitis, but indicate that enterovirus was not detected in the specimen.



Troubleshooting



Factors That Negatively Affect Results

- Improper specimen collection
 - The viral load in the specimen is below the detection limit of the test
 - Performance with other specimen types has not been assessed
- Improper transport or storage of collected specimen
 - Storage and transport conditions are specimen specific
 - Refer to the Package Insert for the appropriate handling instructions
- Improper testing procedure
 - Modification to the testing procedures may alter the performance of the test
 - Careful compliance with the package insert is necessary to avoid erroneous results

INVALID


Test Result **INVALID**

Presence or absence of EV cannot be determined.

SPC/IC does not meet the acceptance criteria, the sample was not properly processed, or PCR was inhibited.

- EV—INVALID
- CIC (SPC/IC)—FAIL
- Probe Check—PASS

Repeat test according to the instructions in the Retest Procedure within the Package Insert.



ERROR

Test Result **ERROR**

Presence or absence of EV cannot be determined.

- EV—NO RESULT
- CIC (SPC/IC)—NO RESULT
- Probe Check—FAIL

Repeat test according to the instructions in the Retest Procedure within the Package Insert.

Test and Analyte Result				Detail	Errors	History
Troubleshoot						
#	Description	Detail		Time		
1	Post-run analysis error	Error 5007 - [EV] probe check failed. Probe check value of 84.19999694824219 for reading number 2 was below the minimum of 98.0		7/6/2007 15:34:50		

Test and Analyte Result				Detail	Errors	History
Assay Name		xpert EV		Version		2
Test Result		ERROR				
Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result		
EV	0.0	0.0	NO RESULT	FAIL		
CIC	0.0	0.0	NO RESULT	PASS		

No Result

Test Result **NO RESULT**

The presence or absence of EV cannot be determined.

- EV—NO RESULT
- CIC (SPC/IC)—NO RESULT
- Probe Check—NA

Possible Causes

A NO RESULT indicates that insufficient data were collected.

- Test was stopped with stop test button
- Electrical failure

Solution

- Secure the power
- Repeat the test with a new cartridge



EV Retest Procedure

1

Discard used cartridge.

Follow your institution's safety guidelines for disposal of cartridges.

2



Obtain a fresh sample.

3



Obtain a new cartridge.

Process the sample per the package insert.

4



Run the test on the system.



Interfering Substances

Positive enterovirus results were obtained even when the highest level of potentially interfering substance was introduced into the assay.

Interfering substance	Concentration	EV C _t
None (Control n = 8)	Not applicable	36.1
Protein (n = 4)	1071 mg / dL	38.2
WBC (n = 4)	7,140 cells / mm ³	37.2
Bloody tap, Specimen 1	2.5% v/v blood	35.9
Bloody tap, Specimen 2	2.5% v/v blood	35.0
Bloody tap, Specimen 3	2.5% v/v blood	35.3
Hemoglobin (n = 4)	3.6 g / dL	36.9

Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
 - Product name
 - Lot number
 - Serial number of the System
 - Error messages (if any)
 - Software version and, if applicable, Computer Service Tag number
- Log your complaint online using the following link <http://www.cepheid.com/us/support> : *Create a Support Case*

Region	Telephone	Technical Support Email
US	+ 1 888 838 3222	techsupport@cepheid.com
Australia and New Zealand	+ 1800 130 821	techsupportANZ@cepheid.com
Brazil and Latin America	+ 55 11 3524 8373	latamsupport@cepheid.com
China	+ 86 021 5406 5387	techsupportchina@cepheid.com
France	+ 33 563 825 319	support@cepheideurope.com
Germany	+ 49 69 710 480 480	support@cepheideurope.com
India, Bangladesh, Bhutan, Nepal, and Sri Lanka	+ 91 11 48353010	techsupportindia@cepheid.com
Italy	+ 39 800 902 567	support@cepheideurope.com
Japan	+ 0120 95 4886	support@japan.cepheid.com
South Africa	+ 27 861 22 76 35	support@cepheideurope.com
United Kingdom	+ 44 3303 332 533	support@cepheideurope.com
Belgium and Netherlands	+33 563 825 3319	support@cepheideurope.com
Other European, Middle East, and African countries	+ 33 563 825 319 + 971 4 253 3218	support@cepheideurope.com



Thank You.



www.Cepheid.com