



Technical Training Xpert[®] Xpress SARS-CoV-2



  In Vitro Diagnostic Medical Device

302-3815 Rev. B November 2022

Training Agenda

- 1 Reagents
- 2 Sample collection
- 3 Kit storage and handling
- 4 Preparing the cartridge
- 5 Quality controls
- 6 Results analysis
- 7 Discussion



Training Objectives

At the end of the training, users will be able to:

- Properly store and handle the Xpert[®] Xpress SARS-CoV-2 kit
- Follow proper laboratory safety precautions
- Collect and transport appropriate specimen
- Prepare a cartridge and run the Xpert[®] Xpress SARS-CoV-2 test
- Report the various software generated results
- Understand the Xpert[®] Xpress SARS-CoV-2 control strategy

The Cepheid Solution



- Detection of SARS-CoV-2
- On-board internal controls for each sample
 - Probe Check Control (PCC)
 - Sample Processing Control (SPC)
- Closed cartridge system minimizes risk of contamination
- On-demand results
- Random access

Intended Use

- The Xpert® Xpress SARS-CoV-2 test is a real-time RT-PCR test intended for the qualitative detection of nucleic acid from the SARS-CoV-2 in nasopharyngeal swab, nasal swab, or nasal wash/aspirate specimen collected from individuals who are suspected of COVID-19 infection.
- Results are for the identification of SARS-CoV-2 RNA. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.
- Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.
- The Xpert® Xpress SARS-CoV-2 test is intended to be performed by trained users in both laboratory and near patient testing settings.

Good Laboratory Practice Review

Personal Protective Equipment (PPE)

- Wear clean lab coats, safety glasses, and gloves
- Change gloves between processing samples

Lab Bench Area

- Clean work surfaces routinely with:
 - ✓ 1:10 dilution of household bleach*
 - ✓ 70% Ethanol Solution
- After cleaning, ensure work surfaces are dry

Specimens, Samples, and Kits Storage

- Store specimens and samples away from kit to prevent contamination

Equipment

- Use filtered pipette tips when recommended
- Follow the manufacturer's requirements for calibration and maintenance of equipment

* Final active chlorine concentration should be 0.5% regardless of the household bleach concentration in your country.

Kit Storage and Handling

Xpert[®] Xpress SARS-CoV-2/Flu/RSV Requirements

GeneXpert[®] Systems

- GeneXpert Dx software **v4.7b** or higher
- Xpertise software **v6.4b** or higher

Test Kits

- XPRSARS-COV2-10

Materials Required but Not Provided

- Nylon flocked swab (Copan P/N 502CS01, 503CS01) or equivalent
- 3mL viral transport medium (Copan P/N 330C) or equivalent
- 3mL 0.9% (w/v) saline
- Nasopharyngeal and Nasal (Copan Part Number 50 units : 305C & 346C) or equivalent

Other Materials

- Personal Protective Equipment (PPE)
- 1:10 dilution bleach
- 70% ethanol or denatured ethanol

Xpert® Xpress SARS-CoV-2/Flu/RSV Kit Contents

Catalog Number	XPRSARS-COV2-10
Cartridges* Per Kit	10
Kit CD	Xpert® Xpress SARS-CoV-2 Assay Definition File (ADF) Xpert® Xpress SARS-CoV-2 Import Instructions Flyer- instructions to access on-line reference materials including the Product Insert
Storage	2-28 °C
Transfer pipettes	10-12



Xpert® Xpress SARS-CoV-2/Flu/RSV Kit Storage and Handling

- Store the Xpert ® Xpress SARS-CoV-2 cartridges and reagents at **2–28°C**
- Follow your institution’s safety procedures for working with chemicals and handling biological samples
- Do not use collection devices that have not been validated by Cepheid
- Open the cartridge lid only when adding the sample, close the lid and proceed with processing
- Start the test within **30 minutes** of adding the sample to the cartridge.

Warnings and Precautions

- Do not shake the cartridge
- Do not use a cartridge... :
 - if it appears wet, has leaked, or if the lid seal appears to have been broken
 - if it appears damaged
 - that has been dropped after removing it from packaging
 - that has been dropped or shaken after you have added the sample
 - that has a damaged reaction tube
 - that has been used; each cartridge is single-use to process one test
 - that is expired
- Do not reuse pipettes

Dispose of cartridges and reagents according to your institution's and country's guidelines for disposal of hazardous materials.

Warnings and Precautions

- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents and require use of standard precautions.
- Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures.
- If national or regional regulations do not provide clear direction on proper disposal, biological specimens and used cartridges should be disposed per WHO [World Health Organization] medical waste handling and disposal guidelines.

Limitations

- Performance characteristics of this test have been established with the specimen types listed in the Intended Use Section only. The performance of this assay with other specimen types or samples has not been evaluated.
- A false negative result may occur if a specimen is improperly collected, transported or handled. False negative results may also occur if inadequate numbers of organisms are present in the specimen.
- As with any molecular test, mutations within the target regions of Xpert Xpress SARS-CoV-2 could affect primer and/or probe binding resulting in failure to detect the presence of virus.
- This test cannot rule out diseases caused by other bacterial or viral pathogens.

For detailed information, refer to the current Instructions For Use

Specimen Collection, Storage and Transport

Specimen Collection

Specimen Type:

nasopharyngeal swab, nasal swab, or nasal wash/ aspirate specimens

Place specimen into 3mL of viral transport medium or 3mL saline



← Nasopharyngeal swab


← Nasal swab

Refer to the CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons Under Investigation (PUIs) for Coronavirus Disease 2019 (COVID-19)

<https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>

Specimen Collection- Nasopharyngeal Swab

1. Insert the swab into either nostril, passing it into the posterior nasopharynx.
2. Rotate swab by firmly brushing against the nasopharynx several times.
3. Remove and place the swab into the tube containing 3mL of viral transport medium or 3mL of saline.
4. Break swab at the indicated break line and cap the specimen collection tube tightly.



Specimen Collection- Nasopharyngeal Swab

Nasopharyngeal Specimen Collection

- 1 Open the package that contains the swab and transport medium tube. Set the tube aside before collecting the specimen.




- 5 Remove the cap from the tube. Insert the swab into the transport medium.




- 2 Open the swab wrapper and remove the swab, taking care not to touch the tip of the swab to any surface.



- 6 Break the swab shaft against the side of the tube at the scoreline. Avoid splashing contents on the skin. Wash with soap and water if exposed.



- 3 Hold the swab in your hand, pinching in the middle of the swab shaft on the scoreline.



- 7 Replace the cap on the tube and close tightly.



- 4 Gently insert the swab into the nostril until you touch the posterior nasopharynx.



Rotate swab several times.

For Xpert Xpress Flu, Xpert Xpress Flu/RSV, and Xpert Xpress SARS-CoV-2/Flu/RSV:

Specimen may be stored for 24 hours at 15-30°C or up to 7 days at 2-8°C.

For Xpert Xpress SARS-CoV-2:

Specimen may be stored for up to 8 hours at 15-30°C or up to 7 days at 2-8°C.

* SWAB/B-100 contains Copan UTM 330C and Copan nylon swab 503CS01



© 2016-2020 Cepheid. All rights reserved.

301-6052, Rev. F October 2020



Specimen Collection- Nasal Swab

1. Insert the nasal swab 1 to 1.5cm into the nostril.
2. Rotate the swab against the inside of the nostril for 3 seconds while applying pressure with a finger to the outside of the nostril.
3. Repeat on the other nostril with the same swab.
4. Remove and place the swab into the tube containing 3mL of viral transport medium or 3mL of saline.
5. Break swab at the indicated break line and cap the specimen collection tube tightly.



Specimen Collection- Nasal Swab

Nasal Swab Specimen Collection


- 1 Open the package that contains the swab and transport medium tube. Set the tube aside before collecting the specimen.



- 2 Open the swab wrapper and remove the swab, taking care not to touch the tip of the swab to any surface.



- 3 Hold the swab in your hand, pinching in the middle of the swab shaft on the scoreline.



- 4 Rotate swab against the inside of the nostril for 3 seconds while applying pressure with a finger to the outside of the nostril.



Do not insert the swabs more than 1-1.5 cm.

- 5 Repeat Step 4 on the other nostril with the same swab.
To avoid specimen contamination, do not touch the swab tip to anything after collecting the specimen.



- 6 Remove the cap from the tube. Insert the swab into the transport medium.



- 7 Break the swab shaft against the side of the tube at the scoreline.
Avoid splashing contents on the skin. Wash with soap and water if exposed.



- 8 Replace the cap on the tube and close tightly.



For Xpert Xpress Flu, Xpert Xpress Flu/RSV, and Xpert Xpress SARS-CoV-2/Flu/RSV:

Specimen may be stored for 24 hours at 15-30°C or up to 7 days at 2-8°C.

For Xpert Xpress SARS-CoV-2:

Specimen may be stored for 8 hours at 15-30°C or up to 7 days at 2-8°C.



Specimen Collection- Nasal Wash/Aspirate

Nasal wash/aspirate specimens can be collected following the user institution standard procedure. Also, refer to the WHO guidelines for the collection of human nasal wash/aspirate specimens.

1. Using a transfer pipette, transfer 600 μ L of the undiluted nasal wash/aspirate specimen into the tube containing 3mL of viral transport medium or 3mL of saline.
2. Cap the tube.

https://www.who.int/influenza/human_animal_interface/virology_laboratories_and_vaccines/guidelines_collection_h5n1_humans/en/

Specimen Transport and Storage

Sample type	Transport and Storage Conditions
Viral Transport Medium or saline containing: nasopharyngeal swab Or nasal swab Or nasal wash/aspirate specimens	 Up to 8 hours  Up to 7 days


Cartridge Preparation



Proper Cartridge Handling Techniques

Correct

- Do not touch the reaction tube
- Keep the cartridge upright after seal has been broken
- Do not tilt when scanning the cartridge



Cartridge Preparation

Xpert® Xpress SARS-CoV-2 Cartridge Preparation

Refer to the package insert for detailed instructions, precautions, and warnings.

For a copy of the SDS, visit www.cepheid.com or www.cepheidinternational.com

Contact information for all Cepheid Technical Support offices is available on our website: www.cepheid.com/en/Customersupport.



- 1 Take one Xpert cartridge for each sample.



- 2 Rapidly invert the tube 5 times.



- 3 Open the cartridge lid.



- 4 Using a clean 300 μ L pipette (supplied), transfer 300 μ L (one draw), of the sample to the opening of the cartridge.



- 5 Close the cartridge lid.



- 6 Start the test within the timeframe specified in the package insert.


Run a Test on GeneXpert[®] Dx

1 Create a test.



Start the test within **30 minutes** after adding the sample to the cartridge.

2 Scan barcode for Patient and/or Sample ID.



Do not click on Manual Entry or Cancel.

3 Scan the cartridge.



Create a Test on GeneXpert Dx Software

4 Complete the fields as required

5 Select the Assay Protocol

Xpert Xpress SARS-CoV-2

6 The module is selected automatically

7 Click on Start Test

8 A green light will flash on the module
Load the cartridge into the module and close the door

The screenshot shows the 'Create Test' window with the following fields and values:


- Patient ID: [Empty]
- Sample ID: [Empty]
- Patient ID 2: [Empty]
- Last Name: [Empty]
- Name: [Empty]
- Select Assay: Xpert Xpress SARS-CoV-2
- Select Module: A3
- Reagent Lot ID*: 16119
- Expiration Date*: 2018/1/17
- Test Type: Specimen
- Sample Type: Other
- Notes: [Empty text area]

Buttons at the bottom: Start Test (highlighted with an orange box and a mouse cursor), Scan Cartridge Barcode.




Run a Test on GeneXpert® Infinity

1 Create a test.




Place the cartridge on the conveyor within **30 minutes** of adding the sample.

2 Scan barcode for Patient and/or Sample ID.



3 Scan the cartridge.



Create a Test on Xpertise Software


4 Complete the fields as required

Order Test - Test Information


Patient ID patientid	
Sample ID sampleid	
Last Name patient	First Name id
Assay Xpert Xpress SARS-CoV-2	
Reagent Lot ID* 12102	Cartridge S/N* 282769448
Expiration Date* 2018/11/04	Priority Normal
Test Type Specimen	
Sample Type Other	Other Sample Type
Notes	

5 Select the Assay Protocol

6 Click on SUBMIT



7 Place the cartridge onto the conveyor belt



Automated Xpert[®] Protocol



Quality Controls

Xpert® Xpress SARS-CoV-2/Flu/RSV Control Strategy CONTROL

- Xpert® Xpress SARS-CoV-2 Quality Controls
 - Each Xpert cartridge is a self-contained test device
 - Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge:
 - Probe Check Controls (PCC)
 - Sample Processing Control (SPC)

Refer to 301-4868 GeneXpert Quality Control Features for All Cepheid Xpert Assays

Refer to 301-4868 GeneXpert® Quality Control Features for all Cepheid Xpert tests.

Internal Quality Controls

- **Probe Check Controls (PCC)**

Before the PCR step, fluorescence signal is measured on all probes and compared with default factory settings to monitor

- reagent rehydration
- probe integrity
- PCR tube filling
- dye stability

- **Sample Processing Controls (SPC)**

SPC ensures that the sample was processed correctly and verifies that sample processing was adequate.

- Verifies adequate extraction and amplification of the sample
- Verifies lysis and detects PCR inhibition
- Must be positive in a negative sample
- Can be positive or negative in a positive sample

Commercially Available External Controls


Vendor	Description	Configuration	Storage
SeraCare AccuPlex™ SARS-CoV-2 Reference Material Kit Catalog # 0505-0126	Positive Control	5 x 1.5mL	2-8°C or -20°C
	Negative Control	5 x 1.5mL	2-8°C or -20°C

1. Open the cartridge lid.
2. Rapidly invert the external control tube 5 times.
3. Using a clean transfer pipette, transfer one draw of the external control sample (300µL) into the large opening (Sample Chamber) in the cartridge.
4. Close cartridge lid.

To minimize degradation of the control material, return any unused sample to the recommended storage conditions immediately after use.

- Many other vendors for quality control material are also available in addition to the one outlined above.
- External controls should be used in accordance with local, state accrediting organizations, as applicable

Result Interpretation



Early Assay Termination

- The Xpert[®] Xpress SARS-CoV-2 test includes an Early Assay Termination (EAT) function which will provide earlier time to results in high titer specimens.
- When SARS-CoV-2 titers are high enough to initiate the EAT function, the SPC amplification curve may not be seen and its results may not be reported.


Results Summary

Result displayed	N2	E	SPC
SARS-CoV-2 POSITIVE	+	+	+/-
	+	-	
SARS-CoV-2 PRESUMPTIVE POS	-	+	+/-
SARS-CoV-2 NEGATIVE	-	-	+
INVALID	-	-	-
ERROR	NO RESULT	NO RESULT	NO RESULT
No Result	NO RESULT	NO RESULT	NO RESULT

SARS-CoV-2 POSITIVE

Test Result

SARS-CoV-2 POSITIVE




The 2019 novel coronavirus (SARS-CoV-2) target nucleic acids are detected.

- The SARS-CoV-2 signal for the N2 nucleic acid target or signals for both nucleic acid targets (N2 and E) have a Ct within the valid range and endpoint above the minimum setting
- SPC: NA; SPC is ignored because coronavirus target amplification occurred
- Probe Check: PASS; all probe check results pass

Test Result

SARS-CoV-2 PRESUMPTIVE POS

SARS-CoV-2 PRESUMPTIVE POS



The 2019 novel coronavirus (SARS-CoV-2) nucleic acids may be present.


Sample should be retested. For samples with a repeated Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.

- The SARS-CoV-2 signal for only the E nucleic acid target has a Ct within the valid range and endpoint above the minimum setting
- SPC: NA; SPC is ignored because a target amplification has occurred.
- Probe Check: PASS; all probe check results pass

SARS-CoV-2 NEGATIVE

Test Result

SARS-CoV-2 NEGATIVE



The 2019 novel coronavirus (SARS-CoV-2) target nucleic acids are not detected.

- The SARS-CoV-2 signals for two nucleic acid targets (N2 and E) do not have a Ct within the valid range and endpoint above the minimum setting
- SPC: PASS; SPC has a Ct within the valid range and endpoint above the minimum setting
- Probe Check: PASS; all probe check results pass

Troubleshooting

Factors That Negatively Affect Results


- Improper specimen collection.
 - The performance of this assay with other specimen types or samples has not been evaluated.
- Inadequate numbers of organisms are present in the specimen.
- Improper transport or storage of collected specimen.
 - Storage and transport conditions are specimen specific.
 - Refer to the Instructions For Use for the appropriate handling instructions.
- Improper testing procedure.
 - Modification to the testing procedures may alter the performance of the test.
 - Careful compliance with the Instructions For Use is necessary to avoid erroneous results.

Reasons to Repeat the Assay

- A **PRESUMPTIVE POS** indicates the 2019 novel coronavirus (SARS-CoV-2) nucleic acids may be present. Only one of the SARS-CoV-2 nucleic acid target was detected (E gene) while the other SARS-CoV-2 nucleic acid target (N2 gene) was not detected.
- An **INVALID** result indicates that the control SPC failed. The sample was not properly processed, PCR is inhibited, or the sample was not properly collected.
- An **ERROR** result could be due to, but not limited to, Probe Check Control failure, system component failure, or the maximum pressure limits were exceeded.
- A **NO RESULT** indicates that insufficient data were collected. For example, cartridge failed integrity test, the operator stopped a test that was in progress, or a power failure occurred.

If an External Control fails to perform as expected, repeat external control test and/or contact Cepheid for assistance.

INVALID Result



SPC does not meet acceptance criteria. Presence or absence of the 2019 novel coronavirus (SARS-CoV-2) nucleic acids cannot be determined.

- SPC: FAIL; SPC and SARS-CoV-2 signals do not have a Ct within valid range and endpoint below minimum setting
- Probe Check – PASS; all probe check results pass

Possible Causes

- Improper sample collection or preparation
- Presence of interfering substances in the sample

Solution

- Repeat the test with a new cartridge



ERROR Result

#	Description	Detail	Time
1	Post-run analysis error	Error 5007: [SCC] probe check failed. Probe check value of 0 for reading number 2 was below the minimum of 33	01/25/15 05:07:22
2	Post-run analysis error	Error 5007: [SPC] probe check failed. Probe check value of 0 for reading number 2 was below the minimum of 222	01/25/15 05:07:22

Presence or absence of the 2019 novel coronavirus (SARS-CoV-2) nucleic acids cannot be determined. Repeat test according to the Retest Procedure in IFU (Section 17.2).

- SARS-CoV-2: NO RESULT
- SPC: NO RESULT
- Probe Check: FAIL; all or one of the probe check results fail

If the probe check passes, the error is caused by the maximum pressure limit exceeding the acceptable range or by a system component failure.

Solution

- Repeat the test with a new cartridge.

NO RESULT

Test Result **NO RESULT**

Test Result	Analyte Result	Detail	Errors	History	Support
Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result	
E	0.0	0	NO RESULT	NA	
N2	0.0	0	NO RESULT	NA	
SPC	0.0	0	NO RESULT	NA	

Presence or absence of the 2019 novel coronavirus (SARS-CoV-2) nucleic acids cannot be determined. A **NO RESULT** indicates that insufficient data were collected. For example, the operator stopped a test that was in progress.

- SARS-CoV-2: NO RESULT
- SPC: NO RESULT
- Probe Check: NA (not applicable)

Possible Causes

A NO RESULT indicates that insufficient data were collected.

- Test was stopped with stop test button
- Electrical failure

Solution

- Secure the power
- Repeat the test with a new cartridge.

Retest Procedure

1

Discard used cartridge

Follow your institution's safety guidelines for disposal of cartridges

2



Obtain the residual specimen, mix according to the Instructions For Use

If the leftover specimen volume is insufficient, or the retest continues to return an INSTRUMENT ERROR or NO RESULT, collect a new specimen.

3




Obtain a new cartridge

Process the specimen per the Instructions For Use

4



Run the test on the System



Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
 - Product name
 - Lot number
 - Serial number of the System
 - Error messages (if any)
 - Software version
- Log your complaint online using the following link
<http://www.cephid.com/en/support>: *Create a Support Case*



Thank You

www.Cepheid.com