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1. Acta Obstet Gynecol Scand. 2008;87(3):260-71. doi: 10.1080/00016340801908759

2. Remington J, et al. Infectious diseases of the fetus and newborn infant. 7th ed. Philadelphia: Elsevier Saunders; 2011.

3. Gonçalves BP, et al. Group B streptococcus infection during pregnancy and infancy: estimates of regional and global burden. Lancet Glob Health. 2022 Apr 28

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1. Ferula A, et al. Bacterial and host determinants of group B streptococcal infection of the neonate and infant. Front Microbiol. 2022 Feb 21;13:820365.

2. Gonçalves BP, et al. Group B streptococcus infection during pregnancy and infancy: estimates of regional and global burden. Lancet Glob Health. 28 de abril de 2022

3. Health P, et al. Neonatal infections: group B streptococcus. BMJ Clin Evid. 2014; 2014: 0323.

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2. Hasperhoven GF, et al. Universal screening versus risk-based protocols for antibiotic prophylaxis during childbirth to prevent early-onset group B streptococcal disease: a systematic review and meta-analysis. BJOG. 2020 May;127(6):680–91.

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3. Verani JR, et al. Division of Bacterial Diseases, National Center for Immunization and Respiratory Diseases, Centers for Disease Control and Prevention (CDC). Prevention of perinatal group B streptococcal disease—revised guidelines from CDC, 2010. MMWR Recomm Rep. 2010 Nov 19;59(RR-10):1–36

4. Wicker E, et al. Group B streptococci: declining incidence in infants in Germany. Pediatr Infect Dis J. 2019 May;38(5):516–9.

5. Angstetra D, et al. Institution of universal screening for Group B streptococcus (GBS) from a risk management protocol results in reduction of early-onset GBS disease in a tertiary obstetric unit. Aust N Z J Obstet Gynaecol. 2007 Oct;47(5):378–82.

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1. Wicker E, et al. Group B streptococci: declining incidence in infants in Germany. *Pediatr Infect Dis J*. 2019 May;38(5):516–9

2. Zimmermann P, et al. Effect of intrapartum antibiotics on the intestinal microbiota of infants: a systematic review. *Arch Dis Child Fetal Neonatal Ed*. 2020 Mar;105(2):201–8

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
1. Di Renzo et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Maternal Fetal Neonatal Med. 2014;1-17. Disponible en: <https://pubmed.ncbi.nlm.nih.gov/25162923>
2. ASM, March 2020, Guidelines for the Detection and Identification of Group B Streptococcus - Revised Guidelines from CDC, 2020

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However, some new developments in the microbiology laboratory and clinical data, highlighted in the new ACOG and ASM guidelines, point out that perhaps the time has arrived for a new update.

Di Renzo et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Maternal Fetal Neonatal Med. 2014:1-17. Disponible en: <https://pubmed.ncbi.nlm.nih.gov/25162923>  
Manuel Rosa-Fraile, Juan-Ignacio Alós. Group B Streptococcus neonatal infections, the ongoing history. 2022.

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1. Helmig R, et al. Diagnostic accuracy of polymerase chain reaction for intrapartum detection of Group B Streptococcus colonization. Acta Obstet Gynecol Scand. 2017 Sep;96(9):1070-1074

2. Young BC, et al. Evaluation of a rapid, real-time intrapartum group B streptococcus assay. Am J Obstet Gynecol. 2011 Oct;205(4):372.e1-6.

3. Zietek M, et al. Intrapartum PCR assay is a fast and efficient screening method for Group B Streptococcus detection in pregnancy. Ginekol Pol. 2020;91(9):549–53

4. Björklund V, et al. Replacing risk-based early-onset-disease prevention with intrapartum group B streptococcus PCR testing. J Matern Fetal Neonatal Med. 2017 Feb;30(3):368-373.

5. Koppes et al. Clinical value of polymerase chain reaction in detecting group B streptococcus during labor. J Obstet Gynaecol Res. 2017 Jun;43(6):996-1000

6. Di Renzo et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Maternal Fetal Neonatal Med. 2014:1-17. Disponible en: <https://pubmed.ncbi.nlm.nih.gov/25162923>

7. ASM, March 2020, Guidelines for the Detection and Identification of Group B Streptococcus - Revised Guidelines from CDC, 2020

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Intrapartum PCR-assay for detection of Group B Streptococci (GBS)

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1. Towers CV, et al. The accuracy of late third-trimester antenatal screening for group B streptococcus in predicting colonization at delivery. Am J Perinatol. 2010 Nov;27(10):785-90.  
2. Björklund V, et al. Replacing risk-based early-onset-disease prevention with intrapartum group B streptococcus PCR testing. J Matern Fetal Neonatal Med. 2017 Feb;30(3):368-373.  
3. Koppes et al. Clinical value of polymerase chain reaction in detecting group B streptococcus during labor. J Obstet Gynaecol Res. 2017 Jun;43(6):996-1000  
4. Helmig RB, et al. Intrapartum PCR-assay for detection of Group B Streptococci (GBS). European Journal of Obstetrics & Gynecology and Reproductive Biology: X. 2019 Oct;4:100081

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**Table 4. Cross-Tabulation of the Antepartum Group B Streptococcus (GBS) Culture Results and the Intrapartum Molecular GBS Test Results**

Antepartum GBS culture result	Intrapartum molecular GBS test result				Total
	Negative	Positive	Invalid <sup>a</sup>	Error <sup>b</sup>	
Negative	664	63	37	54	818
Positive	40	65	4	6	115
Missing	21	10	1	3	35
Total	725	138	42	63	968

<sup>a</sup> Unavailable molecular GBS test results because of polymerase chain reaction inhibition.

<sup>b</sup> Unavailable molecular GBS test results because of significant presence of mucus or manipulation errors in loading the cartridge at the beginning of the study.

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
1. El Helali et al. Diagnostic Accuracy of a Rapid Real-Time Polymerase Chain Reaction Assay for Universal Intrapartum Group B Streptococcus Screening. Clin Infect Dis. (2009) 49 (3): 417-423.

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
PAI

Effect of intrapartum antibiotics on the intestinal microbiota of infants: a systematic review

Petra Zimmermann<sup>1, 2, 3, 4</sup>, Nigel Curtis<sup>2, 3, 4</sup>



Correspondence to Dr Petra Zimmermann, Department of Paediatrics, Faculty of Science and Medicine, University of Fribourg, CH-1700 Fribourg, Switzerland; [petra.zimmermann@mcricri.edu.au](mailto:petra.zimmermann@mcricri.edu.au)

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
1. Lu Zhuang et al. Intestinal Microbiota in Early Life and Its Implications on Childhood Health, Genomics, Proteomics & Bioinformatics, 2019, 17(1), 13–25

2. L. Trasande et al. Infant antibiotic exposures and early-life boy mass. Int. J. Obes. (Lond), 2013

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**Figure 3.** Hierarchical clustering based on composition, at family level, of samples collected at the different times from term infants not exposed to antibiotics and the 4 groups of preterm infants classified as a function of the maternal and/or infant antibiotic administration. Every sample group is associated with its own aggregate representation at family level.

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1. Arboleya et al. Intestinal Microbiota Development in Preterm Neonates and Effect of Perinatal Antibiotics, J of Pediatr 2015

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Immediate GBS results you can trust,  
at the right time — as early as 30 minutes\*

**Coverage, plus Accuracy, plus Peace of mind**  
That's the PCRplus advantage. From Cepheid.

### Impact on Patient Pathway

**Challenge**  
Increased Threat to Newborns Due to Inaccurate Antepartum or Risk-Based Screening

- Antepartum Negative**  
High Risk of GBS Status Change Following Screening at 35-37 Weeks Pregnancy<sup>1</sup>
- Antepartum Positive**  
Incorrect Administration of IAP<sup>2</sup> May Result in EOD<sup>3</sup> and May Impact the Correct Assembly of Fetaloid Flora<sup>4</sup>
- GBS Status Unknown**  
Some Pregnant Women with Unknown GBS Status During Labour<sup>5</sup>

→ Poor outcomes for newborn baby

**Solution**  
Fast & Accurate Intrapartum PCR Testing for GBS as early as 30 minutes\*

- Correct Identification of GBS Status During Delivery
- Evidence-Based IAP & Patient Management
- Reduced Infection Risks & Overall Hospital Costs<sup>6</sup>

### Clinical Performance

Xpert Xpress GBS test results and estimated performance by specimen collection type.

Specimen Collection Type	Sensitivity (95% CI)	Specificity (95% CI)
Intrapartum vaginal/rectal	93.5% (85.7–97.2)	95.5% (93.9–96.7)
Antepartum vaginal/rectal	88.1% (81.1–92.8)	95.6% (93.5–97.0)

\* For positive GBS results with Early Assay Turnaround (EAT), 43 minutes for negative GBS results.  
<sup>1</sup> IAP: Intrapartum Antibiotic Prophylaxis, EOD: Early-Onset Disease, AMR: Antimicrobial Resistance

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IAP. Such test used intrapartum, to be clinically useful, should gather several crucial characteristics as (a) a short turnaround time, (b) accuracy with high sensitivity and specificity, not inferior to 90–95% and 95–98% respectively,<sup>1</sup>

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IAP. Such test used intrapartum, to be clinically useful, should gather several crucial characteristics as (a) a short turnaround time, (b) accuracy with high sensitivity and specificity, not inferior to 90–95% and 95–98% respectively, (c) easiness to perform and to interpret by labor and delivery staff with a minimum of skill and training, and (d) availability at all times 24 hours a day, seven days a week.<sup>1</sup>

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1. Davis HD, et al. Multicenter Study of a Rapid Molecular-Based Assay for the Diagnosis of Group B Streptococcus Colonization in Pregnant Women. C. Infectious Disease. 2004; 30: 1129-35.  
2. Melin P. Neonatal group B streptococcal disease: from pathogenesis to preventive strategies. Clin Microbiol Infect 2011 Sep;17(9):1294-303.  
3. Van Dyke MK, et al. Evaluation of Universal Antenatal Screening for Group B Streptococcus. N Engl J Med 2009;360:2626-36.  
4. Centers for Disease Control and Prevention (CDC). Disparities in universal prenatal screening for group B streptococcus--North Carolina, 2002-2003. MMWR Morb Mortal Wkly Rep. 2005 Jul 22;54(28):700-3.

5. Saari et al. Antibiotic exposure in infancy and risk of being overweight in the first 24 months of life. Pediatrics 2015; 135: 617-626  
6. Di Renzo GC, et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Matern Fetal Neonatal Med. 2015 May;28(7):766-82.  
7. Picchiassi E, et al. Intrapartum test for detection of Group B Streptococcus colonization during labor. J Matern Fetal Neonatal Med. 2017 Aug 31:1-8.

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[www.cepheid.com](http://www.cepheid.com)