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1. Acta Obstet Gynecol Scand. 2008;87(3):260-71. doi: 10.1080/00016340801908759

2. Remington J, et al. Infectious diseases of the fetus and newborn infant. 7th ed. Philadelphia: Elsevier Saunders; 2011.

3. Gonçalves BP, et al. Group B streptococcus infection during pregnancy and infancy: estimates of regional and global burden. Lancet Glob Health. 2022 Apr 28

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1. Ferula A, et al. Bacterial and host determinants of group B streptococcal infection of the neonate and infant. Front Microbiol. 2022 Feb 21;13:820365.

2. Gonçalves BP, et al. Group B streptococcus infection during pregnancy and infancy: estimates of regional and global burden. Lancet Glob Health. 28 de abril de 2022

3. Health P, et al. Neonatal infections: group B streptococcus. BMJ Clin Evid. 2014; 2014: 0323.

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1. Prevention of Early-onset Neonatal Group B Streptococcal Disease: Green-top Guideline No. 36. BJOG. 2017 Nov;124(12):e280–305 [último acceso 15 de junio de 2022]
2. Hasperhoven GF, et al. Universal screening versus risk-based protocols for antibiotic prophylaxis during childbirth to prevent early-onset group B streptococcal disease: a systematic review and meta-analysis. BJOG. 2020 May;127(6):680–91.

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3. Verani JR, et al. Division of Bacterial Diseases, National Center for Immunization and Respiratory Diseases, Centers for Disease Control and Prevention (CDC). Prevention of perinatal group B streptococcal disease—revised guidelines from CDC, 2010. MMWR Recomm Rep. 2010 Nov 19;59(RR-10):1–36
4. Wicker E, et al. Group B streptococci: declining incidence in infants in Germany. Pediatr Infect Dis J. 2019 May;38(5):516–9.
5. Angstetra D, et al. Institution of universal screening for Group B streptococcus (GBS) from a risk management protocol results in reduction of early-onset GBS disease in a tertiary obstetric unit. Aust N Z J Obstet Gynaecol. 2007 Oct;47(5):378–82.

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2. Hasperhoven GF, et al. Universal screening versus risk-based protocols for antibiotic prophylaxis during childbirth to prevent early-onset group B streptococcal disease: a systematic review and meta-analysis. BJOG. 2020 May;127(6):680–91.

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
1. Verani JR, et al. Division of Bacterial Diseases, National Center for Immunization and Respiratory Diseases, Centers for Disease Control and Prevention (CDC). Prevention of perinatal group B streptococcal disease—revised guidelines from CDC, 2010. MMWR Recomm Rep. 2010 Nov 19;59(RR-10):1–36

2. Wicker E, et al. Group B streptococci: declining incidence in infants in Germany. Pediatr Infect Dis J. 2019 May;38(5):516–9.

3. Angstetra D, et al. Institution of universal screening for Group B streptococcus (GBS) from a risk management protocol results in reduction of early-onset GBS disease in a tertiary obstetric unit. Aust N Z J Obstet Gynaecol. 2007 Oct;47(5):378–82.

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
1. Wicker E, et al. Group B streptococci: declining incidence in infants in Germany. *Pediatr Infect Dis J.* 2019 May;38(5):516–9

2. Zimmermann P, et al. Effect of intrapartum antibiotics on the intestinal microbiota of infants: a systematic review. *Arch Dis Child Fetal Neonatal Ed.* 2020 Mar;105(2):201–8

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
1. Di Renzo et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Maternal Fetal Neonatal Med. 2014;1-17. Disponible à l'adresse : <https://pubmed.ncbi.nlm.nih.gov/25162923>
2. ASM, March 2020, Guidelines for the Detection and Identification of Group B Streptococcus - Revised Guidelines from CDC, 2020

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However, some new developments in the microbiology laboratory and clinical data, highlighted in the new ACOG and ASM guidelines, point out that perhaps the time has arrived for a new update.

Di Renzo et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Maternal Fetal Neonatal Med. 2014:1-17.

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Manuel Rosa-Fraile, Juan-Ignacio Alós. Group B Streptococcus neonatal infections, the ongoing history. 2022.

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1. elmig R, et al. Diagnostic accuracy of polymerase chain reaction for intrapartum detection of Group B Streptococcus colonization. Acta Obstet Gynecol Scand. 2017 Sep;96(9):1070-1074

2. Young BC, et al. Evaluation of a rapid, real-time intrapartum group B streptococcus assay. Am J Obstet Gynecol. 2011 Oct;205(4):372.e1-6.

3. Zietek M, et al. Intrapartum PCR assay is a fast and efficient screening method for Group B Streptococcus detection in pregnancy. Ginekol Pol. 2020;91(9):549-53

4. Björklund V, et al. Replacing risk-based early-onset-disease prevention with intrapartum group B streptococcus PCR testing. J Matern Fetal Neonatal Med. 2017 Feb;30(3):368-373.

5. Koppes et al. Clinical value of polymerase chain reaction in detecting group B streptococcus during labor. J Obstet Gynaecol Res. 2017 Jun;43(6):996-1000

6. Di Renzo et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Maternal Fetal Neonatal Med. 2014:1-17. Disponible à l'adresse : <https://pubmed.ncbi.nlm.nih.gov/25162923>

7. ASM, March 2020, Guidelines for the Detection and Identification of Group B Streptococcus - Revised Guidelines from CDC, 2020

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Intrapartum PCR-assay for detection of Group B Streptococci (GBS)

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1. Towers CV, et al. The accuracy of late third-trimester antenatal screening for group B streptococcus in predicting colonization at delivery. Am J Perinatol. 2010 Nov;27(10):785-90.
2. Björklund V, et al. Replacing risk-based early-onset-disease prevention with intrapartum group B streptococcus PCR testing. J Matern Fetal Neonatal Med. 2017 Feb;30(3):368-373.
3. Koppes et al. Clinical value of polymerase chain reaction in detecting group B streptococcus during labor. J Obstet Gynaecol Res. 2017 Jun;43(6):996-1000
4. Helmig RB, et al. Intrapartum PCR-assay for detection of Group B Streptococci (GBS). European Journal of Obstetrics & Gynecology and Reproductive Biology: X. 2019 Oct;4:100081

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Table 4. Cross-Tabulation of the Antepartum Group B Streptococcus (GBS) Culture Results and the Intrapartum Molecular GBS Test Results

Antepartum GBS culture result	Intrapartum molecular GBS test result				Total
	Negative	Positive	Invalid ^a	Error ^b	
Negative	664	63	37	54	818
Positive	40	65	4	6	115
Missing	21	10	1	3	35
Total	725	138	42	63	968

^a Unavailable molecular GBS test results because of polymerase chain reaction inhibition.

^b Unavailable molecular GBS test results because of significant presence of mucus or manipulation errors in loading the cartridge at the beginning of the study.

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
1. El Helali et al. Diagnostic Accuracy of a Rapid Real-Time Polymerase Chain Reaction Assay for Universal Intrapartum Group B Streptococcus Screening. Clin Infect Dis. (2009) 49 (3): 417-423.

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
Effect of intrapartum antibiotics on the intestinal microbiota of infants: a systematic review

Petra Zimmermann^{1, 2, 3, 4}, Nigel Curtis^{2, 3, 4}



Correspondence to Dr Petra Zimmermann, Department of Paediatrics, Faculty of Science and Medicine, University of Fribourg, CH-1700 Fribourg, Switzerland; petra.zimmermann@mcrl.edu.au

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1. Lu Zhuang et al. Intestinal Microbiota in Early Life and Its Implications on Childhood Health, Genomics, Proteomics & Bioinformatics, 2019, 17(1), 13–25

2. L. Trasande et al. Infant antibiotic exposures and early-life boy mass. Int. J. Obes. (Lond), 2013

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
Figure 3. Hierarchical clustering based on composition, at family level, of samples collected at the different times from term infants not exposed to antibiotics and the 4 groups of preterm infants classified as a function of the maternal and/or infant antibiotic administration. Every sample group is associated with its own aggregate representation at family level.

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1. Arboleya et al. Intestinal Microbiota Development in Preterm Neonates and Effect of Perinatal Antibiotics, J of Pediatr 2015

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Immediate GBS results you can trust,
at the right time — as early as 30 minutes*

Xpert® Xpress GBS

Coverage, plus Accuracy, plus Peace of mind
That's the PCRplus advantage. From Cepheid.

Impact on Patient Pathway

Challenge
Increased Threat to Newborns Due to Inaccurate Antepartum or Risk-Based Screening

Pregnant Woman

Healthcare Professional

Antepartum Negative

Antepartum Positive

GBS Status Unknown

High Risk of GBS Status Change Following Screening at 35-37 Weeks Pregnancy¹

Incorrect Administration of IAP² May Result in EOD³ and May Impact the Correct Assembly of Fetaloid Flora⁴

Some Pregnant Women with Unknown GBS Status During Labour⁵

→ Poor outcomes for newborn baby

Solution
Fast & Accurate Intrapartum PCR Testing for GBS as early as 30 minutes*

Pregnant Woman

Healthcare Professional

Fast PCR Across Lab & Decentralised Settings, 24/7

Correct Identification of GBS Status During Delivery

Evidence-Based IAP² & Patient Management

Reduced Infection Risks & Overall Hospital Costs⁶

Clinical Performance

Xpert Xpress GBS test results and estimated performance by specimen collection type.

Specimen Collection Type	Sensitivity (95% CI)	Specificity (95% CI)
Intrapartum vaginal/rectal	93.5% (85.7–97.2)	95.5% (93.9–96.7)
Antepartum vaginal/rectal	88.1% (81.1–92.8)	95.6% (93.5–97.0)

* For positive GBS results with Early Access Transcription (EAT), 45 minutes for negative GBS results.
¹ IAP: Intrapartum Antibiotic Prophylaxis, EOD: Early-Onset Disease, AMR: Antimicrobial Resistance

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IAP. Such test used intrapartum, to be clinically useful, should gather several crucial characteristics as (a) a short turnaround time, (b) accuracy with high sensitivity and specificity, not inferior to 90–95% and 95–98% respectively,¹

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IAP. Such test used intrapartum, to be clinically useful, should gather several crucial characteristics as (a) a short turnaround time, (b) accuracy with high sensitivity and specificity, not inferior to 90–95% and 95–98% respectively, (c) easiness to perform and to interpret by labor and delivery staff with a minimum of skill and training, and (d) availability at all times 24 hours a day, seven days a week.¹

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1. Davis HD, et al. Multicenter Study of a Rapid Molecular-Based Assay for the Diagnosis of Group B Streptococcus Colonization in Pregnant Women. C. Infectious Disease. 2004; 30: 1129-35.
2. Melin P. Neonatal group B streptococcal disease: from pathogenesis to preventive strategies. Clin Microbiol Infect 2011 Sep;17(9):1294-303.
3. Van Dyke MK, et al. Evaluation of Universal Antenatal Screening for Group B Streptococcus. N Engl J Med 2009;360:2626-36.
4. Centers for Disease Control and Prevention (CDC). Disparities in universal prenatal screening for group B streptococcus--North Carolina, 2002-2003. MMWR Morb Mortal Wkly Rep. 2005 Jul 22;54(28):700-3.

5. Saari et al. Antibiotic exposure in infancy and risk of being overweight in the first 24 months of life. Pediatrics 2015; 135: 617-626
6. Di Renzo GC, et al. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Matern Fetal Neonatal Med. 2015 May;28(7):766-82.
7. Picchiasse E, et al. Intrapartum test for detection of Group B Streptococcus colonization during labor. J Matern Fetal Neonatal Med. 2017 Aug 31:1-8.

Questions/ réponses

Merci

www.cepheid.com