

Xpert® Xpress SARS-CoV-2

Rapid, accurate detection of SARS-CoV-2



The Need

- SARS-CoV-2 infections are associated with increased **morbidity, mortality, cost** and pose an ever-present **threat to public health** globally.^{1,2}
- Prior emergence of SARS-CoV-2 **variants caused new surges in cases and death.**²
- Testing required to identify the breakthrough infections and inform patient care, especially for the **high-risk population.**³

The Impact

- Support clinicians with on-demand, timely, and accurate results.
- Significant reduction in time-to-results reduces isolation time and total costs while optimizing capacity.⁴
- Enable patients to receive timely and appropriate treatment.

The Solution

The Xpert **Xpress** SARS-CoV-2 test is a real-time PCR test intended for the qualitative detection of nucleic acid from the SARS-CoV-2 in nasopharyngeal swab, nasal swab, or specimen collected from individuals who are suspected of COVID-19 infection.

The Xpert **Xpress** SARS-CoV-2 test provides:

- Fast, accurate results in as early as 30 minutes.*
- Two gene targets for the detection of SARS-CoV-2.
- Rapid sample-to-answer testing with actionable results from a single sample.

¹ WHO Coronavirus (COVID-19) Dashboard –Last accessed 22nd Sep 2023

² <https://news.un.org/en/story/2023/05/1136367>

³ Peeling RW, Heymann DL, Teo YY, Garcia PJ. Diagnostics for COVID-19: moving from pandemic response to control. *Lancet*. 2022 Feb 19;399(10326):757-768. doi: 10.1016/S0140-6736(21)02346-1. Epub 2021 Dec 20. PMID: 34942102; PMCID: PMC8687671

⁴ Fistera D, Kikull K, Risse J, Herrmann A, Brachmann M, Kill C. Point-of-care PCR testing of SARS-CoV-2 in the emergency department: Influence on workflow and efficiency. *PLoS One*. 2023 Aug 3;18(8) <https://pubmed.ncbi.nlm.nih.gov/37535577>

* With early assay termination for positives, otherwise, the full test runtime is 45 minutes.



Xpert® Xpress SARS-CoV-2

Product Reference Sheet — CE-IVD

Test Reagent Kit	Xpert Xpress SARS-CoV-2	
Catalog Number	XPRSARS-COV2-10	
Technology	Real-time RT-PCR	
Targets	N2 – nucleocapsid gene E – envelope protein gene	
Batch or On-Demand	On-Demand	
Minimum Batch Size	1	
Sample Types	Specimen Collection: Nasopharyngeal swab, nasal swab Transport Media: UTM/VTM or Saline	
Sample Extraction	Automated/integrated	
Precision Pipetting	Not required	
Turnaround Time	As soon as 30 minutes for positives* and approximately 45 minutes for negatives	
Hands-on Time	< 1 minute	
Controls: Process	Sample Processing Control	
Controls: Probe Function/Detection	Probe Check Control	
	Positive Percent Agreement	Negative Percent Agreement
Clinical Evaluation	97.8% (95% CI: 88.4%–99.6%)	95.6% (95% CI: 85.2%–98.8%)
	<i>Testing performed with 45 positives and 45 negatives</i>	
Sample Storage	15–30 °C for up to 8 hours or 2–8 °C for up to 7 days until testing is performed	
Kit Storage	2–28 °C	
Commercial Controls	Refer to Xpert Xpress SARS-CoV-2 Package Insert or contact Cepheid Technical Support	

* With Early Assay Termination (EAT) for positive results.
Refer to current package insert for complete details.

CE-IVD. *In Vitro* Diagnostic Medical Device. May not be available in all countries.

CORPORATE HEADQUARTERS

904 Caribbean Drive
Sunnyvale, CA 94089 USA

TOLL FREE +1.888.336.2743
PHONE +1.408.541.4191
FAX +1.408.541.4192

EUROPEAN HEADQUARTERS

Vira Solelh
81470 Maurens-Scopont France

PHONE +33.563.82.53.00
FAX +33.563.82.53.01
EMAIL cepheid@cepheideurope.fr

www.Cepheidinternational.com

© 2023–2024 Cepheid. 3344-02

