



Technical Training Xpert® Xpress SARS-CoV-2/Flu/RSV

For Use with GeneXpert® Dx or
GeneXpert® Infinity Systems

*Catalog Number (Ex: GXGBSLBXC-10)
For CE-IVD Only*



302-5322 Rev. C November 2022



Training Agenda

- 1 Reagents
- 2 Sample collection
- 3 Kit storage and handling
- 4 Preparing the cartridge
- 5 Quality controls
- 6 Results analysis
- 7 Discussion



Training Objectives

At the end of the training, users will be able to:

- Properly store and handle the Xpert[®] Xpress SARS-CoV-2/Flu/RSV kit
- Follow proper laboratory safety precautions
- Collect and transport appropriate specimen
- Prepare a cartridge and run the Xpert[®] Xpress SARS-CoV-2/Flu/RSV test
- Report the various software generated results
- Understand the Xpert[®] Xpress SARS-CoV-2/Flu/RSV control strategy

The Cepheid Solution



- Detection of SARS-CoV-2, FluA, FluB, RSV RNA
- On-board internal controls for each sample
 - Probe Check Control (PCC)
 - Sample Processing Control (SPC)
- Closed cartridge system minimizes risk of contamination
- EAT (Early Assay Termination for SARS-CoV-2 ADF only)
- On-demand results
- Random access

Intended Use

- The Xpert® Xpress SARS-CoV-2/Flu/RSV test is a multiplexed real-time RT-PCR test intended for the simultaneous, qualitative detection and differentiation of SARS-CoV-2, influenza A, influenza B, and respiratory syncytial virus (RSV) viral RNA in either nasopharyngeal swab, nasal swab or nasal wash/aspirate specimens collected from individuals suspected of respiratory viral infection. Clinical signs and symptoms of respiratory viral infection due to SARS-CoV-2, influenza, and RSV can be similar.
- Results are for the simultaneous detection and differentiation of SARS-CoV-2, influenza A virus, influenza B virus and RSV RNA in clinical specimens. Positive results are indicative of the presence of the identified virus, but do not rule out bacterial infection or co-infection with other pathogens not detected by the test.
- Negative results do not preclude SARS-CoV-2, influenza A virus, influenza B virus and/or RSV infection and should not be used as the sole basis for treatment or other patient management decisions. Negative results must be combined with clinical observations, patient history, and/or epidemiological information

Good Laboratory Practice Review

Personal Protective Equipment (PPE)

- Wear clean lab coats, safety glasses, and gloves
- Change gloves between processing samples

Lab Bench Area

- Clean work surfaces routinely with:
 - ✓ 1:10 dilution of household bleach*
 - ✓ 70% Ethanol Solution
- After cleaning, ensure work surfaces are dry

Specimens, Samples, and Kits Storage

- Store specimens and samples away from kit to prevent contamination

Equipment

- Use filtered pipette tips when recommended
- Follow the manufacturer's requirements for calibration and maintenance of equipment

* Final active chlorine concentration should be 0.5% regardless of the household bleach concentration in your country.

Specimen Collection, Storage and Transport

Specimen Collection

Specimen Type:

nasopharyngeal swab, nasal swab, or nasal wash/ aspirate specimens

Place specimen into 3mL of viral transport medium or 3mL saline



← Nasopharyngeal swab


← Nasal swab

Refer to the CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons Under Investigation (PUIs) for Coronavirus Disease 2019 (COVID-19)

<https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>

Specimen Collection- Nasopharyngeal Swab

1. Insert the swab into either nostril, passing it into the posterior nasopharynx.
2. Rotate swab by firmly brushing against the nasopharynx several times.
3. Remove and place the swab into the tube containing 3mL of viral transport medium or 3mL of saline.
4. Break swab at the indicated break line and cap the specimen collection tube tightly.



Specimen Collection- Nasopharyngeal Swab

Nasopharyngeal Specimen Collection


- 1 Open the package that contains the swab and transport medium tube. Set the tube aside before collecting the specimen.



- 2 Open the swab wrapper and remove the swab, taking care not to touch the tip of the swab to any surface.



- 3 Hold the swab in your hand, pinching in the middle of the swab shaft on the scoreline.




- 4 Gently insert the swab into the nostril until you touch the posterior nasopharynx. Rotate swab several times.



- 5 Remove the cap from the tube. Insert the swab into the transport medium.



- 6 Break the swab shaft against the side of the tube at the scoreline. Avoid splashing contents on the skin. Wash with soap and water if exposed.



- 7 Replace the cap on the tube and close tightly.



For Xpert Xpress Flu, Xpert Xpress Flu/RSV, and Xpert Xpress SARS-CoV-2/Flu/RSV:

Specimen may be stored for 24 hours at 15-30°C or up to 7 days at 2-8°C.

For Xpert Xpress SARS-CoV-2:

Specimen may be stored for up to 8 hours at 15-30°C or up to 7 days at 2-8°C.

* SWAB/B-100 contains Copan UTM 330C and Copan nylon swab 503CS01


© 2016-2020 Cepheid. All rights reserved.

301-6052, Rev. F October 2020



Specimen Collection- Nasal Swab

1. Insert the nasal swab 1 to 1.5cm into the nostril.
2. Rotate the swab against the inside of the nostril for 3 seconds while applying pressure with a finger to the outside of the nostril.
3. Repeat on the other nostril with the same swab.
4. Remove and place the swab into the tube containing 3mL of viral transport medium or 3mL of saline.
5. Break swab at the indicated break line and cap the specimen collection tube tightly.



Specimen Collection- Nasal Swab

Nasal Swab Specimen Collection


- 1 Open the package that contains the swab and transport medium tube. Set the tube aside before collecting the specimen.



- 2 Open the swab wrapper and remove the swab, taking care not to touch the tip of the swab to any surface.



- 3 Hold the swab in your hand, pinching in the middle of the swab shaft on the scoreline.



- 4 Rotate swab against the inside of the nostril for 3 seconds while applying pressure with a finger to the outside of the nostril.



Do not insert the swabs more than 1-1.5 cm.

- 5 Repeat Step 4 on the other nostril with the same swab.
To avoid specimen contamination, do not touch the swab tip to anything after collecting the specimen.



- 6 Remove the cap from the tube. Insert the swab into the transport medium.



- 7 Break the swab shaft against the side of the tube at the scoreline.
Avoid splashing contents on the skin. Wash with soap and water if exposed.



- 8 Replace the cap on the tube and close tightly.



For Xpert Xpress Flu, Xpert Xpress Flu/RSV, and Xpert Xpress SARS-CoV-2/Flu/RSV:

Specimen may be stored for 24 hours at 15-30°C or up to 7 days at 2-8°C.

For Xpert Xpress SARS-CoV-2:

Specimen may be stored for 8 hours at 15-30°C or up to 7 days at 2-8°C.




Specimen Collection- Nasal Wash/Aspirate

- Nasal wash/aspirate specimens can be collected following the user institution's standard procedure. Also, refer to the WHO guidelines for the collection of human nasal wash/aspirate specimens.

https://www.who.int/influenza/human_animal_interface/virology_laboratories_and_vaccines/guidelines_collection_h5n1_humans/en/

- Using a clean transfer pipette, transfer 600 μL of the sample into the tube containing 3 mL of viral transport medium or 3 mL of saline and then cap the tube.

Specimen Transport and Storage

Sample type	Transport and Storage Conditions
<p>3mL Viral Transport Medium containing nasopharyngeal swab, nasal swab, or nasal wash/aspirate</p>	 ≤ 24 hours
<p>3mL saline containing nasopharyngeal swab, nasal swab, or nasal wash/aspirate</p>	 ≤ 48 hours
<p>3mL Viral Transport Medium or 3mL saline containing nasopharyngeal swab, nasal swab, or nasal wash/aspirate</p>	 ≤ 7 days

Kit Storage and Handling

Xpert® Xpress SARS-CoV-2/Flu/RSV Requirements

GeneXpert® Systems

- GeneXpert Dx software **v4.7b** or higher
- Xpertise software **v6.4b** or higher

Test Kits

- XPCOV2/FLU/RSV-10

Materials Required but Not Provided

- Nylon flocked swab (Copan P/N 502CS01, 503CS01) or equivalent
- 3mL viral transport medium (Copan P/N 330C) or equivalent
- 3mL 0.9% (w/v) saline
- Nasopharyngeal and Nasal (Copan Part Number 50 units : 305C & 346C) or equivalent

Other Materials

- Personal Protective Equipment (PPE)
- 1:10 dilution bleach
- 70% ethanol or denatured ethanol


Xpert® Xpress SARS-CoV-2/Flu/RSV Kit Contents

Catalog Number	XPCOV2/FLU/RSV-10
Cartridges* Per Kit	10
Storage	2-28 °C
Flyer	Directions to locate the Instructions For Use and Quick Reference Instructions on www.cepheid.com www.cepheid.com/coronavirus-resources
Transfer pipettes	10-12

The kit also includes printed copies of the Quick Reference Instructions, which should only be used with the GeneXpert® Xpress System.



Cartridges contain chemically hazardous substances-please see Instructions for Use and Safety Data Sheet for more detailed information.



Xpert[®] Xpress SARS-CoV-2/Flu/RSV Kit Storage and Handling

- Store test kits at 2-28°C. Do not use expired cartridges.
- Each single-use cartridge is used to process one test. Do not reuse processed cartridges.
- Do not open a cartridge until ready to use.
- Start the test within 30 minutes of adding the sample to the cartridge.
- To avoid cross contamination during sample handling steps, change gloves between samples.

Cartridge Preparation




Warnings and Precautions

- Do not shake the cartridge
- Do not use a cartridge... :
 - if it appears wet, has leaked, or if the lid seal appears to have been broken
 - if it appears damaged
 - that has been dropped after removing it from packaging
 - that has been dropped or shaken after you have added the sample
 - that has a damaged reaction tube
 - that has been used; each cartridge is single-use to process one test
 - that has expired
- Do not reuse pipettes
- Do not reuse swabs

Proper Cartridge Handling Techniques

Correct

- Do not touch the reaction tube
- Keep the cartridge upright after seal has been broken
- Do not tilt when scanning the cartridge



Xpert® Xpress SARS-CoV-2/Flu/RSV

Cartridge Preparation

Sample Qualification – check if all items below are present:

1. Transport media containing swab (if applicable)
2. Patient name or identifier on the tube
3. Cartridges and transport media are within the expiration date

Good Laboratory Practices:

- Wear clean gloves, lab coats, and safety glasses.
- Change gloves between samples.
- Clean work surface with 1:10 dilution of bleach followed by 70% ethanol solution.

Xpert® Cartridge Preparation

- Xpert Xpress SARS-CoV-2
- Xpert Xpress SARS-CoV-2/Flu/RSV

Refer to the package insert for detailed instructions, precautions, and warnings.

For a copy of the SDS, visit www.cepheid.com or www.cepheidinternational.com

Contact information for all Cepheid Technical Support offices is available on our website: www.cepheid.com/en/CustomerSupport.



1 Take one Xpert cartridge for each sample.



2 Rapidly invert the tube 5 times.



3 Open the cartridge lid.



4 Using a clean 300 µL pipette (supplied), transfer 300 µL (one draw) of the sample to the opening of the cartridge.



5 Close the cartridge lid.




6 Start the test within the timeframe specified in the package insert.


Xpert[®] Xpress SARS-CoV-2/Flu/RSV Cartridge Preparation



1
Take one Xpert cartridge for each sample.



2
Rapidly invert the tube 5 times.



3
Open the cartridge lid.



4
Using a clean 300 μ L pipette (supplied), transfer 300 μ L (one draw) of the sample to the cartridge.



5
Close the cartridge lid.

6

Start the test within the timeframe specified in the Instructions For Use.


Run a Test on GeneXpert[®] Dx

1 Create a test.



Start the test within **30 minutes** after adding the sample to the cartridge.

2 Scan barcode for Patient and/or Sample ID.



Do not click on Manual Entry or Cancel.

3 Scan the cartridge.



Create a Test on GeneXpert[®] Dx Software

4 Complete the fields as required

5 Select the Assay Protocol

6 The module is selected automatically

7 Click on Start Test

8 A green light will flash on the module
Load the cartridge into the module and close the door

Create Test

Patient ID
Sample ID
Patient ID 2
Last Name

Select Assay
Select Module A3
Reagent Lot ID* 16119 Expiration Date* 2018/1/17
Test Type Specimen
Sample Type Other
Notes


Start Test Scan Cartridge Barcode

Assay	Count
<None>	
<None>	
Xpert Xpress_SARS-CoV-2_Flu_RSV	1
Xpert Xpress_SARS-CoV-2_Flu	1
Xpert Xpress_SARS-CoV-2	1




Run a Test on GeneXpert® Infinity

1 Create a test.



Place the cartridge on the conveyor within **30 minutes** of adding the sample.

2 Scan barcode for Patient and/or Sample ID.



3 Scan the cartridge.



Create a Test on Xpertise Software

4 Complete the fields as required

Order Test - Test Information

Patient ID	patientid
Sample ID	sampleid
Last Name	patient
First Name	id

5 Select the Assay Protocol

Assay Protocol

Xpert Xpress SARS-CoV-2	
-------------------------	--

Reagent Lot ID*
12102

Expiration Date*
2018/11/04

<None>	
Xpert Xpress_SARS-CoV-2_Flu_RSV	1
Xpert Xpress_SARS-CoV-2_Flu	1
Xpert Xpress_SARS-CoV-2	1

Test Type
Specimen

Sample Type
Other


Other Sample Type

Notes

6 Click on SUBMIT



7 Place the cartridge onto the conveyor belt



Waste Disposal

Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents and require use of standard precautions.

Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures.

If national or regional regulations do not provide clear direction on proper disposal, biological specimens and used cartridges should be disposed per WHO [World Health Organization] medical waste handling and disposal guidelines.

Quality Controls

Xpert[®] Xpress SARS-CoV-2/Flu/RSV Control Strategy CONTROL

- Xpert[®] Xpress SARS-CoV-2/Flu/RSV Quality Controls
 - Each Xpert cartridge is a self-contained test device
 - Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge:
 - Probe Check Controls (PCC)
 - Sample Processing Control (SPC)

Refer to 301-4868 GeneXpert[®] Quality Control Features for all Cepheid Xpert tests.

Xpert® Xpress SARS-CoV-2/Flu/RSV Cartridge Controls

- **Probe Check Controls (PCC)**

Before the PCR step, fluorescence signal is measured on all probes and compared with default factory settings to monitor

- reagent rehydration
- probe integrity
- PCR tube filling
- dye stability

- **Sample Processing Controls (SPC)**

SPC ensures that the sample was processed correctly and verifies that sample processing was adequate.

- Verifies adequate extraction and amplification of the sample
- Verifies lysis and detects PCR inhibition
- Must be positive in a negative sample to be a valid test
- Can be positive or negative in a positive sample

Commercially Available External Controls

Zeptomatrix	Description	Configuration	Storage
NATFRC-6C	NATtrol Flu/RSV/SARS-CoV-2 Positive Control	6 x 0.5mL	2-8°C or -20°C
NATCV9-6C	Coxsackievirus A9 Negative Control	6 x 0.5mL	2-8°C or -20°C

1. Open the cartridge lid.
2. Rapidly invert the external control tube 5 times.
3. Using a clean transfer pipette, transfer one draw (300µl) of the external control sample into the large opening (Sample Chamber) in the cartridge.
4. Close cartridge lid.

To minimize degradation of the control material, return any unused sample to the recommended storage conditions immediately after use.

- Many other vendors for quality control material are also available in addition to the one outlined above.
- External controls should be used in accordance with local, state accrediting organizations, as applicable

Result Interpretation

Assay Targets

- SARS-CoV-2
- Flu A1
- Flu A2
- Flu B
- RSV
- SPC

Early Assay Termination

- The Xpert® Xpress SARS-CoV-2 test mode includes an Early Assay Termination (EAT) function which will provide earlier time to results in high titer specimens if the signal from the SARS-CoV-2 target reaches a predetermined threshold before the full 45 PCR cycles have been completed.
- When SARS-CoV-2 titers are high enough to initiate the EAT function, the SPC amplification curve may not be seen and its results may not be reported.

Results Summary SARS-CoV-2 ADF

Result displayed	SARS-CoV-2	SPC
SARS-CoV-2 POSITIVE	+	+/-
SARS-CoV-2, NEGATIVE	-	+
INVALID	-	-
ERROR	NO RESULT	NO RESULT
No Result	NO RESULT	NO RESULT

Results Summary SARS-CoV-2 and Flu ADF

Result displayed	SARS-CoV-2	Flu A1	Flu A2	Flu B	SPC
Influenza A POSITIVE	-	+	+/-	-	+/-
Influenza A POSITIVE	-	+/-	+	-	+/-
Influenza B POSITIVE	-	-	-	+	+/-
SARS-CoV-2 POSITIVE	+	-	-	-	+/-
SARS-CoV-2, Flu A, Flu B, NEGATIVE	-	-	-	-	+
INVALID	-	-	-	-	-
ERROR	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT
No Result	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT


Results Summary SARS-CoV-2, Flu, and RSV ADF

Result displayed	SARS - CoV-2	Flu A1	Flu A2	Flu B	RSV	SPC
Influenza A POSITIVE	-	+	+/-	-	-	+/-
Influenza A POSITIVE	-	+/-	+	-	-	+/-
Influenza B POSITIVE	-	-	-	+	-	+/-
RSV POSITIVE	-	-	-	-	+	+/-
SARS-CoV-2 POSITIVE	+	-	-	-	-	+/-
SARS-CoV-2, Flu A, Flu B, RSV NEGATIVE	-	-	-	-	-	+
INVALID	-	-	-	-	-	-
ERROR	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT
No Result	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT	NO RESULT

SARS-CoV-2 POSITIVE

Test Result **SARS-CoV-2 POSITIVE**

Test Result	Analyte Result	Detail	Errors	History	Messages	Support
Analyte Name	Ct	EndPt	Interpretation Result	Reason	Analyte Result	Probe Check Result
SARS-CoV-2	30.9	57	POS		POS	PASS
SPC	30.2	14	PASS		NA	PASS




- SARS-CoV2 target RNA are detected
- SPC: NA; SPC is ignored because target amplification occurred
- Probe Check: PASS; all probe check results pass

SARS CoV2 Negative, Influenza A Positive, Influenza B Negative

Test Result	SARS-CoV-2 NEGATIVE; Flu A POSITIVE; Flu B NEGATIVE
-------------	---


Test Result	Analyte Result	Detail	Errors	History	Messages	Support
Analyte Name	Ct	EndPt	Interpretation Result	Reason	Analyte Result	Probe Check Result
SARS-CoV-2	0.0	2	NEG	No Ct	NEG	PASS
Flu A 1	30.3	576	POS		POS	PASS
Flu A 2	34.1	272	POS		POS	PASS
Flu B	0.0	7	NEG	No Ct	NEG	PASS
SPC	29.2	164	PASS		NA	PASS

- SARS-CoV-2 not detected
Flu A target RNA detected;
Flu B target RNA not detected;
- SPC is ignored because target amplification occurred
- Probe Check: PASS; all probe check results pass



SARS CoV2 Negative, Influenza A Negative, Influenza B Positive

Test Result	Analyte Result	Detail	Errors	History	Messages	Support
Analyte Name	Ct	EndPt	Interpretation Result	Reason	Analyte Result	Probe Check Result
SARS-CoV-2	0.0	-6	NEG	No Ct	NEG	PASS
Flu A 1	0.0	-2	NEG	No Ct	NEG	PASS
Flu A 2	0.0	-1	NEG	No Ct	NEG	PASS
Flu B	28.5	1059	POS		POS	PASS
SPC	29.0	181	PASS		NA	PASS




- SARS-CoV-2 target RNA not detected;
Flu A target RNA not detected;
Flu B target RNA detected;
- SPC is ignored because target amplification occurred
- Probe Check: PASS; all probe check results pass

Test Result	SARS-CoV-2 NEGATIVE; Flu A NEGATIVE; Flu B NEGATIVE; RSV POSITIVE
-------------	--

SARS CoV2 Negative, Influenza A Negative, Influenza B Negative, RSV Positive


Test Result	Analyte Result	Detail	Errors	History	Messages	Support
Analyte Name	Ct	EndPt	Interpretation Result	Reason	Analyte Result	Probe Check Result
SARS-CoV-2	0.0	2	NEG	No Ct	NEG	PASS
Flu A 1	0.0	1	NEG	No Ct	NEG	PASS
Flu A 2	0.0	3	NEG	No Ct	NEG	PASS
Flu B	0.0	2	NEG	No Ct	NEG	PASS
RSV	29.0	389	POS		POS	PASS
SPC	29.6	157	PASS		NA	PASS



- SARS-CoV-2 not detected
Flu A target RNA not detected
Flu B target RNA not detected
RSV target RNA detected
- SPC is ignored because target amplification occurred
- Probe Check: PASS; all probe check results pass

Test Result SARS-CoV-2 NEGATIVE;
Flu A NEGATIVE;
Flu B NEGATIVE;
RSV NEGATIVE

SARS CoV2 Negative, Influenza A Negative, Influenza B Negative, RSV Negative



- SARS-CoV-2 not detected
Flu A target RNA not detected
Flu B target RNA not detected
RSV target RNA not detected
- SPC: PASS; SPC has a Ct within the valid range and endpoint above the minimum setting
- Probe Check: PASS; all probe check results pass

Limitations

- Performance of the Xpert Xpress SARS-CoV-2/Flu/RSV test has only been established in nasopharyngeal and nasal swab specimens. Use of the Xpert Xpress SARS-CoV-2/Flu/RSV test with other specimen types has not been assessed and performance characteristics are unknown.
- Nasal wash/aspirate specimens are considered acceptable specimen types for use with the Xpert Xpress SARS-CoV-2/Flu/RSV test but performance with these specimen types has not been established .
- As with any molecular test, mutations within the target regions of the Xpert Xpress SARS-CoV-2/Flu/RSV test could affect primer and/or probe binding resulting in failure to detect the presence of virus or the virus being detected less predictably.
- This test cannot rule out diseases caused by other bacterial or viral pathogens.
- The performance of this test was validated using the procedures provided in this package insert only. Modifications to these procedures may alter the performance of the test.
- Erroneous test results might occur from improper specimen collection; failure to follow the recommended sample collection, handling, and storage procedures; technical error; or sample mix-up. Careful compliance with the instructions in this insert is necessary to avoid erroneous results.
- False negative results may occur if virus is present at levels below the analytical limit of detection.

Limitations (continued)

- Negative results do not preclude SARS-CoV-2, influenza or RSV infection and should not be used as the sole basis for treatment or other patient management decisions.
- Results from the Xpert Xpress SARS-CoV-2/Flu/RSV test should be correlated with the clinical history, epidemiological data, and other data available to the clinician evaluating the patient.
- Viral nucleic acid may persist *in vivo*, independent of virus viability. Detection of analyte target(s) does not imply that the corresponding virus(es) are infectious or are the causative agents for clinical symptoms.
- This test has been evaluated for use with human specimen material only.
- This test is a qualitative test and does not provide the quantitative value of detected organism present.
- This test has not been evaluated for monitoring treatment of infection.
- This test has not been evaluated for screening of blood or blood products for the presence of SARS-CoV-2, influenza or RSV.
- The effect of interfering substances has only been evaluated for those listed within the labeling. Interference by substances other than those described can lead to erroneous results.

Limitations (continued)

- Results from analytical studies with contrived co-infected samples showed potential for competitive interference when SARS-CoV-2, influenza or RSV was present at 1X LoD levels.
- Cross-reactivity with respiratory tract organisms other than those described herein can lead to erroneous results.
- Recent patient exposure to FluMist® or other live attenuated influenza vaccines may cause inaccurate positive results.
- As the Xpert Xpress SARS-CoV-2/Flu/RSV test does not differentiate between the N2 and E gene targets, the presence of other coronaviruses in the B lineage, *Betacoronavirus* genus, including SARS-CoV-1 may cause a false positive result. None of these other coronaviruses is known to currently circulate in the human population.
- This test is not intended to differentiate RSV subgroups, influenza A subtypes or influenza B lineages. If differentiation of specific RSV or influenza subtypes and strains is needed, additional testing, in consultation with state or local public health departments, is required.
- Specimen transport media that contain guanidine thiocyanate (GTC) may interfere with the test causing false negative results.

Troubleshooting


Factors That Negatively Affect Results

- Improper specimen collection.
 - The performance of this assay with other specimen types or samples has not been evaluated.
- Inadequate numbers of organisms are present in the specimen.
- Improper transport or storage of collected specimen.
 - Storage and transport conditions are specimen specific.
 - Refer to the Instructions For Use for the appropriate handling instructions.
- Improper testing procedure.
 - Modification to the testing procedures may alter the performance of the test.
 - Careful compliance with the Instructions For Use is necessary to avoid erroneous results.

Test Result	INVALID
-------------	---------

INVALID Result

Test Result	Analyte Result	Detail	Errors	History	Messages	Support
Analyte Name	Ct	EndPt	Interpretation Result	Reason	Analyte Result	Probe Check Result
SARS-CoV-2	0.0	1	NEG	No Ct	INVALID	PASS
Flu A 1	0.0	0	NEG	No Ct	INVALID	PASS
Flu A 2	0.0	1	NEG	No Ct	INVALID	PASS
Flu B	0.0	-1	NEG	No Ct	INVALID	PASS
RSV	0.0	0	NEG	No Ct	INVALID	PASS
SPC	0.0	3	FAIL	No Ct	FAIL	PASS



SPC does not meet acceptance criteria. Presence or absence of the target RNA cannot be determined.

- SPC: FAIL;
- SARS-CoV-2, Flu A, Flu B, RSV signals do not have a Ct within valid range and endpoint below minimum setting
- Probe Check – PASS; all probe check results pass

Possible Causes

- Improper sample collection or preparation
- Presence of interfering substances in the sample

Solution

- Repeat the test with a new cartridge

Test Result	ERROR
-------------	--------------

ERROR Result

Test Result				Analyte Result				Detail				Errors				History				Messages				Support			
Troubleshoot																											
#	Description	Detail										Time															
1	Operation terminated	Error 2125: Termination Error - Insufficient Volume: 17, 50, 0, 18										07/13/20 17:21:51															

Presence or absence of the target RNA cannot be determined.

- SARS-CoV-2: NO RESULT
- Flu A: NO RESULT
- Flu B: NO RESULT
- RSV: NO RESULT
- SPC: NO RESULT
- Probe Check: FAIL; all or one of the probe check results fail

If the probe check passes, the error is caused by the maximum pressure limit exceeding the acceptable range or by a system component failure.

Solution

- Repeat the test with a new cartridge.

NO RESULT

Test Result	Analyte Result	Detail	Errors	History	Messages	Support
Analyte Name	Ct	EndPt	Interpretation Result	Reason	Analyte Result	Probe Check Result
SARS-CoV-2	0.0	0	NO RESULT		NO RESULT	NA
Flu A 1	0.0	0	NO RESULT		NO RESULT	NA
Flu A 2	0.0	0	NO RESULT		NO RESULT	NA
Flu B	0.0	0	NO RESULT		NO RESULT	NA
RSV	0.0	0	NO RESULT		NO RESULT	NA
SPC	0.0	0	NO RESULT		NO RESULT	NA

Presence or absence of the target RNA cannot be determined.

A **NO RESULT** indicates that insufficient data was collected. For example, the operator stopped a test that was in progress.

Possible Causes

- A NO RESULT indicates that insufficient data was collected.
 - Test was stopped with stop test button
 - Electrical failure

Solution

- Secure the power
- Repeat the test with a new cartridge.

Reasons to Repeat the Assay

- An **INVALID** result indicates that the control SPC failed. The sample was not properly processed, PCR is inhibited, or the sample was not properly collected.
- An **ERROR** result could be due to, but not limited to, Probe Check Control failure, system component failure, no sample added, or the maximum pressure limits were exceeded.
- A **NO RESULT** indicates that insufficient data were collected. For example, cartridge failed integrity test, the operator stopped a test that was in progress, or a power failure occurred.

Retest Procedure

1

Discard used cartridge.
Follow your institution's
safety guidelines for
disposal of cartridges.


2



If enough specimen is
available, re-test from
original specimen collection
tube.

If insufficient specimen is
available, a new specimen
must be collected.

3




Obtain a new cartridge.

Process the specimen per
the Instructions For Use.

4



Run the test on the system.



Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
 - Product name
 - Lot number
 - Serial number of the System
 - Error messages (if any)
 - Software version
- Log your complaint online using the following link
<http://www.cephid.com/en/support>: *Create a Support Case*



Thank You

www.Cepheid.com