

# Xpert® Bladder Cancer Detection Xpert® Bladder Cancer Monitor

*Two non-invasive urine tests that deliver accurate and actionable results*



## The Need

Hematuria is one of the most common symptoms of bladder cancer.<sup>1</sup> Almost all patients with bladder cancer present hematuria, but less than 10% of patients with hematuria are diagnosed with bladder cancer.<sup>1</sup>

Patients diagnosed with non-muscle invasive bladder cancer (NMIBC) have high recurrence and progression, requiring frequent follow-ups.<sup>11</sup>

The European Association of Urology (EAU) and the American Urological Association (AUA) recommend cystoscopy and histology to diagnose patients with hematuria and cystoscopy for surveillance of patients diagnosed with NMIBC and not treated with radical cystectomy.<sup>2,3</sup> The standard methods are invasive, painful, relatively expensive, operator-dependent, and often imperfect.<sup>3</sup>

## The Impact

- Fast and noninvasive test to detect or exclude bladder cancer and determine the need for further diagnostic follow-up.<sup>4,5</sup>
- Enables use of clinically relevant and readily available urine specimens to avoid unnecessary procedures and discomfort.<sup>4,6</sup>
- Helps detect early onset and recurrence of bladder cancer before it is detectable by cystoscopy.<sup>5,7,8</sup>

## The Solution

The **Xpert® Bladder Cancer Detection** is intended as an aid to detect the presence of bladder cancer in patients with hematuria suspected of having bladder cancer.<sup>9</sup>

The **Xpert® Bladder Cancer Monitor** is intended as an aid to monitor the recurrence of bladder cancer in patients previously diagnosed with bladder cancer.<sup>10</sup>

The Xpert® Bladder Cancer Detection and Xpert Bladder Cancer Monitor tests are non-invasive, easy to implement, fast, and painless biomarker tests with high negative predictive values.<sup>9,10</sup> Both tests utilize voided urine specimens to measure the level of five target mRNAs (ABL1, CRH, IGF2, UPK1B, ANXA10) in a self-contained cartridge and can be performed in any size urology lab without a traditional PCR lab setup.<sup>9,10</sup>

- 1 Tae Jeong Oh, Ji Yong Lee, Yangyei Seo, Min A. Woo, Jae Sung Lim, Yong Gil Na, Ki Hak Song, Bo-Ram Bang, Justin Junguek Lee, Ju Hyun Shin, Sungwhan An, Evaluation of Sensitive Urine DNA-Based PENK Methylation Test for Detecting Bladder Cancer in Patients with Hematuria, *The Journal of Molecular Diagnostics*, Volume 25, Issue 9, 2023, Pages 646-654, ISSN 1525-1578, <https://doi.org/10.1016/j.jmoldx.2023.05.003>.
- 2 EAU Guidelines. Edn. presented at the EAU Annual Congress Milan 2023.
- 3 Chang SS, Boorjian SA, Chou R. Diagnosis and treatment of non-muscle invasive bladder cancer: AUA/SUO guideline. *AUA/SUO Joint Guideline*: Published 2016; Amended 2020 [Internet]. <https://www.auanet.org/>. 2020 [cited 2021 Jun 21]. Available from: <https://www.auanet.org/guidelines/guidelines/bladder-cancer-non-muscle-invasive-guideline>
- 4 Palou J, Brausi M, Catto JWF. Management of Patients with Normal Cystoscopy but Positive Cytology or Urine Markers. *Eur Urol Oncol*. 2020;3(4):548–54.
- 5 Valenberg FJPV, Hiar AM, Wallace E, Bridge JA, Mayne DJ, Beqaj S, Sexton WJ, Lotan Y, Weizer AZ, Jansz GK, Stenzl A, Danella JF, Cline KJ, Williams MB, Montgomery S, David RD, Harris R, Klein EW, Bradford TJ, Wolk FN, Westenfelder KR, Trainer AF, Richardson TA, Egerdie RB, Goldfarb B, Zadra JA, Lu X, Simon IM, Campbell SA, Bates MP, Higuchi RG, Witjes JA. Validation of an mRNA-based Urine Test for the Detection of Bladder Cancer in Patients with Haematuria. *Eur Urol Oncol*. 2021 Feb;4(1):93-101. doi: 10.1016/j.euo.2020.09.001. Epub 2020 Sep 28. PMID: 33004290.
- 6 Warrick JI, Sjö Dahl G, Kaag M, Raman JD, Merrill S, Shuman L, et al. Intratumoral heterogeneity of bladder cancer by molecular subtypes and histologic variants. *Eur Urol*. 2019 Jan;75(1):18–22.
- 7 Cowan B, Klein E, Jansz K, Westenfelder K, Bradford T, Peterson C, et al. Longitudinal follow-up and performance validation of an mRNA-based urine test (Xpert® Bladder Cancer Monitor) for surveillance in patients with non-muscle-invasive bladder cancer. *BJU Int*. 2021 Dec;128(6):713–21.
- 8 Soorojebally Y, Neuzillet Y, Roumiguié M, Lamy PJ, Allory Y, Descotes F, Ferlicot S, Kassab-Chahmi D, Oudard S, Rébillard X, Roy C, Leuret T, Rouprêt M, Audenet F. Urinary biomarkers for bladder cancer diagnosis and NMIBC follow-up: a systematic review. *World J Urol*. 2023 Feb;41(2):345-359. doi: 10.1007/s00345-022-04253-3. Epub 2023 Jan 2. PMID: 36592175.
- 9 Xpert Bladder Cancer Detection Instructions for Use. 301-2414.
- 10 Xpert Bladder Cancer Monitor Instructions for Use. 301-5933.
- 11 Witjes JA. Follow-up in non-muscle invasive bladder cancer: facts and future. *World J Urol*. 2021 Nov;39(11):4047-4053. doi: 10.1007/s00345-020-03569-2. Epub 2020 Dec 26. PMID: 33367941; PMCID: PMC8571151.



# Xpert® Bladder Cancer Detection

## Product Reference Sheet — CE-IVD

<b>Test Reagent Kit</b>	Xpert Bladder Cancer Detection
<b>Catalog Number</b>	<b>CE-IVD</b> GXBLAD-CD-CE-10
<b>Technology</b>	Real-time RT-PCR
<b>Targets</b>	Five mRNA targets UPK1B, IGF2, CRH, ANXA10, ABL1
<b>Batch or On-Demand</b>	On-demand
<b>Minimum Batch Size</b>	1
<b>Sample Type</b>	Voided Urine (not first morning void)
<b>Sample Volume</b>	4 ml
<b>Sample Extraction</b>	Automated/integrated
<b>Precision Pipetting</b>	Not Required
<b>Off-board Sample Preparation Time</b>	Approximately 5 minutes
<b>TAT</b>	Approximately 90 minutes
<b>Controls: Sample Adequacy Control</b>	ABL1
<b>Controls: Probe Function/Detection</b>	Probe Check Control (PCC)
<b>Controls: Internal Control</b>	Armored RNA®

### Xpert Bladder Cancer Detection Performance vs. Cystoscopy/Histology

<i>(N=895 study subjects)</i>	<b>Overall</b>	<b>Low Grade</b>	<b>High Grade</b>
<b>Sensitivity</b>	75.8%	52.2%	88.4%
<b>Specificity</b>	84.6%	N/A	N/A
<b>Negative Predictive Value</b>	97.8%		
<b>Positive Predictive Value</b>	28.1%		
<b>System &amp; Software</b>	<b>GeneXpert Dx System</b> GeneXpert Dx software version 4.7b or higher	<b>GeneXpert Infinity</b> Xpertise software version 6.4b or higher	
<b>Sample Stability</b>	Urine samples in Xpert Urine Transport Reagent tubes are stable up to 7 days at 2–28 °C		
<b>Kit Storage</b>	2–28 °C		
<b>Additional Required Materials</b>	Xpert Urine Transport Reagent Kit	<b>Catalog Number:</b> GXUTR-CE-30	
<b>Commercial Controls</b>	Refer to Instructions for Use (IFU) or Contact Cepheid Technical Support		



# Xpert® Bladder Cancer Monitor

## Product Reference Sheet — CE-IVD

<b>Test Reagent Kit</b>	Xpert Bladder Cancer Monitor
<b>Catalog Number</b>	<b>CE-IVD</b> GXBLAD-CM-CE-10
<b>Technology</b>	Real-time RT-PCR
<b>Targets</b>	Five mRNA targets UPK1B, IGF2, CRH, ANXA10, ABL1
<b>Batch or On-Demand</b>	On-demand
<b>Minimum Batch Size</b>	1
<b>Sample Type</b>	Voided Urine (not first morning void)
<b>Sample Volume</b>	4 ml
<b>Sample Extraction</b>	Automated/integrated
<b>Precision Pipetting</b>	Not Required
<b>Off-board Sample Preparation Time</b>	Approximately 5 minutes
<b>TAT</b>	Approximately 90 minutes
<b>Controls: Sample Adequacy Control</b>	ABL1
<b>Controls: Probe Function/Detection</b>	Probe Check Control (PCC)
<b>Controls: Internal Control</b>	Armored RNA®

### Xpert Bladder Cancer Monitor Performance vs. Cystoscopy/Histology

(N=255 study subjects)	<b>Overall</b>	<b>Low Grade</b>	<b>High Grade</b>
<b>Sensitivity</b>	75.0%	63.2%	84.0%
<b>Specificity</b>	80.6%	N/A	N/A
<b>Negative Predictive Value</b>	93.9%		
<b>Positive Predictive Value</b>	44.6%		
<b>System &amp; Software</b>	<b>GeneXpert Dx System</b> GeneXpert Dx software version 4.7b or higher	<b>GeneXpert Infinity</b> Xpertise software version 6.4b or higher	
<b>Sample Stability</b>	Urine samples in Xpert Urine Transport Reagent tubes are stable up to 7 days at 2–28 °C		
<b>Kit Storage</b>	2–28 °C		
<b>Additional Required Materials</b>	Xpert Urine Transport Reagent Kit	<b>Catalog Number:</b> GXUTR-CE-30	
<b>Commercial Controls</b>	Refer to Instructions for Use (IFU) or Contact Cepheid Technical Support		

CE-IVD. In Vitro Diagnostic Medical Device. May not be available in all countries. Not available in the United States.

#### CORPORATE HEADQUARTERS

904 Caribbean Drive  
Sunnyvale, CA 94089 USA

TOLL FREE +1.888.336.2743  
PHONE +1.408.541.4191  
FAX +1.408.541.4192

#### EUROPEAN HEADQUARTERS

Vira Soleth  
81470 Maurens-Scopont France

PHONE +33.563.82.53.00  
FAX +33.563.82.53.01  
EMAIL [cepheid@cepheideurope.fr](mailto:cepheid@cepheideurope.fr)

[www.Cepheidinternational.com](http://www.Cepheidinternational.com)

© 2024 Cepheid. 3354-02

