



Xpert® Xpress SARS-CoV-2/Flu/RSV

Actionable Rapid Respiratory Results to Meet the Challenges of the Ongoing Pandemic



The Need

- Detection and differentiation of SARS-CoV-2, flu, and RSV which all present with similar symptoms
- Insights for effective triage and patient management
- Risk assessment for patients and staff while resuming traditional healthcare procedures
- Identification of potential cases of co-infections during respiratory season
- Current flu and SARS-CoV-2 therapies need definitive detection for each pathogen



The Solution

- Actionable detection of SARS-CoV-2, Flu A, Flu B, and RSV in as little as 25 minutes*
- Actionable on-demand results with one sample collection
- Optimization of GeneXpert® system module capacity by combining two **Xpress** test cartridges into one
- Standardization of results between the central lab and near patient testing sites



The Impact

- Empower patients with rapid results to help keep their loved ones and the community safe
- Assist hospitals to effectively manage limited resources
- Provide labs with single-test confidence to accelerate workflows
- Support clinicians with on-demand, timely, and accurate results

Coverage, plus Accuracy, plus Peace of mind

That's the **PCRplus** advantage. From Cepheid.

Simple Implementation into Existing Workflows

As Simple as 1, 2, 3 — Minimal Hands-on Time

1

Obtain specimen stored in appropriate transport medium



2

Transfer sample to cartridge



3

Scan and insert cartridge



Xpert® Xpress SARS-CoV-2/Flu/RSV is a qualitative real-time PCR test that provides accurate detection and differentiation of SARS-CoV-2, Flu A, Flu B, and RSV in as little as 25 minutes.*

Catalog Information

Xpert Xpress SARS-CoV-2/Flu/RSV	10 tests	XPCOV2/FLU/RSV-10
--	----------	-------------------

COLLECTION KITS#

Cepheid Nasopharyngeal Collection Kits	100 units	SWAB/M-100
	100 units	SWAB/B-100
Cepheid Nasal Collection Kits	100 units	SWAB/F-100

SAMPLE TYPES

Specimen Collection	Nasopharyngeal swabs, nasal swabs or nasal wash/aspirates
Transport Media [^]	UTM/VTM or saline

EXTERNAL CONTROLS#

AVAILABLE FROM ZEPTOMETRIX

ZeptoMetrix NATrol Flu/RSV/SARS-CoV-2 Positive Control Kit	6 vials	NATFRC-6C
ZeptoMetrix NATrol Flu/RSV/SARS-CoV-2 Negative Control Kit	6 vials	NATCV9-6C

CE-IVD. *In Vitro* Diagnostic Medical Device. Not available in all countries.

* For positives when running SARS-CoV-2 only on Xpert **Xpress** SARS-CoV-2/Flu/RSV; otherwise, the test will have a runtime of approximately 36 minutes.

[^] Specimen transport media that contain guanidine thiocyanate (GTC) may interfere with the test causing false negative results.

Recommended but not required materials, please see package insert 302-5159 (for GeneXpert Dx and GeneXpert Infinity Systems) for complete details.

CORPORATE HEADQUARTERS

904 Caribbean Drive
Sunnyvale, CA 94089 USA

TOLL FREE +1.888.336.2743
PHONE +1.408.541.4191
FAX +1.408.541.4192

EUROPEAN HEADQUARTERS

Vira Solelh
81470 Maurens-Scopont France

PHONE +33.563.82.53.00
FAX +33.563.82.53.01
EMAIL cepheid@cepheideurope.fr

www.Cepheidinternational.com

© 2022–2023 Cepheid. 3194-03