

Xpert Check v2.0 Software Release Notes

This document contains new information, known limitations, and clarifications not described in the *Xpert Check Package Insert*.

New Features

Xpert Check software version 2.0 includes the following new features or improvements:

- **Decoupling Firmware Version Binding**— Xpert Check now allows a range of firmware versions per module model type. Also added validation for modules based on module model number and firmware version.
- **Update to Database Security** - You can now run Xpert Check with Product software using the latest SQL passwords as per Windows guidelines to connect to the Database. Older SQL passwords are still supported, and the passwords are now encrypted for added security.
- **Update to Contact Information**—The Contact Information info such as serial number and Country will now be validated for accuracy and consistency. Added dropdown menu for Country field.
- **Updated to Reporting**— Reports have also been updated to show “N.A.” for module tests that have been skipped instead of blank fields. The table layout has been updated to accommodate larger module serial numbers.
- **Update to Java 11** – Xpert Check updated to Java 11 to bring in the latest security features, and performance benefits.

Known Issues

No known limitations that may impact overall system performance.