



Upgrade today.

More answers,
More targets,
More insights,

Get **more** than results.

GeneXpert® System

With 10-Color Technology:

Hybrid Module Configuration


10-color modules enable your GeneXpert system to run advanced versions of tests, future multiplexed tests and syndromic panels.

GeneXpert II, IV, XVI Systems with both 6-color and 10-color modules*

10-color modules are calibrated with 4 additional dyes for detection of 40+ identified targets.


New tests, new technology.
No need to buy a new instrument, just upgrade!

Only Dx software 6.5 or higher is validated to run hybrid GeneXpert systems.



Your Cepheid representative will support you in customizing your GeneXpert fleet to your laboratory's needs and keep delivering high-quality testing to your patients.


It has never been easier to get more answers...



1 Daisy Chain

To expand your GeneXpert system at full capacity, add a 10-color instrument


Ensure readiness for Xpert® MTB/XDR and future Xpert tests while keeping a small footprint and saving on extra computer, LIS connection, or validation.



2 Expand Capacity

To maximize your GeneXpert system's remaining capacity, add 10-color modules

Keep your testing capacity aligned with your needs and adopt the newest 10-color technology to ensure the readiness of your fleet for Xpert MTB/XDR or future Xpert tests.



3 Convert Modules

To upgrade a 6-color system that's at capacity, swap its 6-color modules for 10-color

Ensure the readiness of your fleet for Xpert MTB/XDR or future Xpert tests. Option expected in 2024.

* 10-color hybrids are not available for GeneXpert systems R1 generation, GeneXpert Infinity systems, GeneXpert Xpress systems or GeneXpert systems running Microsoft Windows 7 operating system.

IVD. In Vitro Diagnostic Medical Device. May not be available in all countries

CORPORATE HEADQUARTERS

904 Caribbean Drive
Sunnyvale, CA 94089 USA

TOLL FREE +1.888.336.2743
PHONE +1.408.541.4191
FAX +1.408.541.4192

EUROPEAN HEADQUARTERS

Vira Solelh
81470 Maurens-Scopont France

PHONE +33.563.82.53.00
FAX +33.563.82.53.01
EMAIL Cepheid@cepheideurope.fr

www.Cepheidinternational.com

© 2023 Cepheid. 5115-01